IKA ULTRA-TURRAX® T 25 digital

Betriebsanleitung	DE	5
Operating instructions	EN	12
Mode d'emploi	FR	19
Instrucciones de manejo	ES	26
Handleiding	NL	33
Istruzioni per l'uso	IT	40
Driftsanvisning	SV	47
Driftsinstruks	DA	54
Driftsvejledning	NO	61
Käyttöohje	FI	68
Instruções de serviço	PT	75
Instrucja obsługi	PL	82
Návod k provozu	CS	89
Használati utasítás	HU	96
Navodilo za delovanje	SL	103
Návod na použitie	SK	110
Kasutusjuhend	ET	117
Lietošanas instrukcija	LV	124
Darbo instrukcija	LT	131

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Contents

	Page
CE-Declaration of conformity	3
Safety instructions	12
Correct use	13
Unpacking	14
Useful information	14
Drive	14
Assembly of boom and stand	15
Working with the unit	15
Maintenance and cleaning	15
Accessories	16
Permitted dispersion tools	16
Error codes	17
Warranty	18
Technical data	18

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - splashing liquids
 - body parts, hair, clothing and jewellery getting caught.
- The voltage stated on the nameplate must correspond to the mains voltage.


Fig. 3: Points of danger

- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts.
- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids etc.
- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- There may be electrostatic discharges between the medium and the dispersing instrument shaft which could pose a direct danger.
- The appliance is not suitable for manual operation.
- Set up the stand in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that the stand does not start to move.
- The agitated vessels used for stirring have to be secured. Consider on a good stability of the entire structure.
- Secure the stirring vessel against twisting.
- Glass vessels must always be secured with a clamp to prevent them spinning. When working with flasks elastic couplers must be used to prevent fracture of the glass.
- Note the operating instructions of the dispersing tool and accessories.
- Only dispersing elements approved by **IKA** may be used.
- Please observe the permitted speed for the dispersing element. Never set higher speeds.
- Do not use the appliance without a dispersing element.
- Use the dispersing tool always inside the stirring vessel.
- Check that the turning handles are secure and tighten if necessary.
- Make certain that the unit is set at the lowest speed before commissioning; otherwise, the unit will begin running at the speed last set. Gradually increase the speed.
- Reduce the speed if the medium splashes out of the vessel be-

cause the speed is too high.

- Before switching on the dispersing instrument make sure that its shaft is immersed min. 45 mm in the medium to prevent the medium from splashing out.
- The distance between the dispersion tool and the vessel bottom should not be less than 10 mm.
- In the event of unbalance or unusual noises, switch off the appliance immediately. Replace the dispersing element. If there is no difference after the change of the dispersing tool, return it to the dealer or the manufacturer along with a description of the fault.
- Do not touch rotating parts during operation.
- Please note that the dispersing element and the journal bearings can become extremely hot during use.
- Do not cover the ventilation slots on the drive in order to ensure adequate cooling of the drive.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do not process any flammable or combustible materials.
- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA** if you have any questions.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- Meter powder not too close to the flange. Powder can be blown away by air turbulences of the drive.
- The appliance does not start up again automatically following a cut in the power supply.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Always disconnect the plug before fitting accessories.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.
- The gasket and bearings are made of PTFE and rustproof steel; the following points should therefore be noted: Chemical reac-

tions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

Correct use

► Application

When used in combination with one of our recommended dispersing elements, the drive unit is a high-speed dispersing and emulsifying unit capable of handling free-flowing and liquid media in batches.

- | | |
|--------------------|--|
| - Production of: | Emulsions
Dispersions
Wet crushing |
| - Operating modes: | on stand |

► Range of use

- | | |
|----------------|-----------|
| - Laboratories | - Schools |
| - Pharmacies | |

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

Unpacking

► Unpack

- Please unpack the device carefully
- In the case of any damage a fact report must be set immediately (post, rail or forwarder)

► Delivery scope

- Disperser **T 25 digital**
- Extension arm
- Hexagon socket screw
- Tool set for dispersion tools serie S25
- Operating instruction

improvements in the obtainable fineness; the energy expended serves merely to increase the temperature of the medium.


Fig. 4

Useful information

Dispersion is the dissolution and diffusion of a solid, liquid or gaseous phase in a liquid that is not consolute with that phase.

The rotor/stator principle

Due to the high rotation speed of the rotor, the medium to be processed is automatically drawn axially into the dispersion head and then forced radially through the slots in the rotor/stator arrangement. The high accelerations acting on the material produce extremely strong shear and thrust forces. In addition, high turbulence occurs in the shear gap between rotor and stator, which provides optimum mixing of the suspension.

The dispersion effectiveness is heavily dependent on the product of the shear gradient and the time the particles spend in the shear zone. The optimum range for the circumferential velocity of the rotor/stator arrangement is 6-24 m/s.

A processing time of a few minutes is usually sufficient to produce the desired fineness. Long processing times bring only insignificant

Drive unit

The drive unit opens up a wide range of possibilities for dispersion technology in the conventional laboratory with a respective performance output of approximately 300 watts at 24,000 1/min.

The speed is set by using the differently adjusting wheel. The motor speed can be read off the LED display. For example, a value of 13.6 corresponds to 13,600 min⁻¹.

The motor speed can be recorded using an **IKA** Tachometer DZM control.o, which can be connected after removing the cover (3) on the boom side of the device (see fig. 1). Data is transferred to the PC using the RS 232 interface on the DZM control.o.

Assembly of boom and stand

The delivered extension arm is mounted as follows (overview fig. 1):

- Insert extension arm (pos. 1) into flange
- Screw in the cylinder head screw (pos. 2)
- Firmly tighten cylinder head screw (pos. 2) using the A/f 4 right angled screwdriver.

Vibrations may loosen the cylinder head screw. Therefore, for your own safety, periodically check that the extension arm is secure. If necessary, retighten the cylinder head screw.

To ensure safe operation the driving unit is fastened to the plate stand R 1826 (pos. 8) by means of a cross sleeve (pos. 7).

In order to enhance the stability of the mechanical arrangement, the driving unit has to be mounted as close as possible to the stand pillar.

Working with the unit

Also pay attention to the ambient conditions listed in the "Technical data".

Assembling the dispersion tool (overview fig. 1)

- Open the knurled screw (pos. 10) so that the thread does not protrude into the hole.
- Insert the dispersing element (9) from the bottom into the drive unit as far as the stop. The shaft will audibly click into place when a small amount of force is applied (sprung ball bearing).
- **Attention:** Now secure the dispersing element by screwing the knurled screw tight (pos. 10).

The filling level of the medium may not be less than about 55 mm.

The entire unit may also be arranged somewhat off centre in order to avoid any undesired air admission caused by the strong rotational turbulence.

The adjustment wheel (pos. 4) is set at the lowest speed before the unit is switched on: value 1 on the adjustment dial corresponds to 3,400 min⁻¹.

By using the adjustment wheel (pos. 4) the speed may be infinitely varied to meet the requirements of the medium to be processed. The idle speed for the drive unit can be adjusted in the range 3,400 and 24,000 min⁻¹; the value will be shown on the LED display.

If these conditions are met the unit is ready for service when the mains plug has been plugged in.

The driving unit is switched on by means of the On/Off- switch (pos. 6).

Maintenance and cleaning

Drive unit: The drive unit is maintenance-free but not immune from wear. The carbon brushes of the motor wear down over time.

For cleaning disconnect the main plug.

To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Dispersion tools: Consider please the manuals associated to the dispersion tools.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see **www.ika.de**, spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Accessories

- DZM control.o Rpm counter
- R 1826 Plate stand
- R 182 Cross sleeve
- RH 3 Tensioner

Permitted dispersion tools

Dispersing tools serie S25...


Abbreviations used for dispersion tools serie S25...

- N : PTFE bearing
- NK : Ball bearing with PTFE bearing
- KR : Ball bearing with rotary shaft seal (suitable for a vacuum up to 50 mbar)
- KV : Ball bearing with axial face seal (suitable for vacuum up to 1 mbar, 6 bar excess pressure)
- D : without seal
- KS : synthetic material

10G, 14G, 18G, 19G, 25G, 25F:

Stator/shaft tube diameter G: thick; F: fine

Speed range : upto 24,000 min⁻¹

Designation	Shaft material	
S25N-8G	Stainless steel	
S25N-10G	Stainless steel	
S25N-10G-VS	Stainless steel	
S25N-18G	Stainless steel	
S25KR-18G	Stainless steel	
S25KV-18G	Stainless steel	
S25NK-19G	Stainless steel	
S25N-25G	Stainless steel	
S25KR-25G	Stainless steel	
S25KV-25G	Stainless steel	
S25N-25F	Stainless steel	
S25KR-25F	Stainless steel	
S25KV-25F	Stainless steel	
S25KV-25G-IL	Stainless steel	
S25KV-25F-IL	Stainless steel	
S25D-10G-KS	Synthetic mat.	
S25D-14G-KS	Synthetic mat.	

For applications and further information, please see the operating instructions of the dispersing element.

Only use the dispersion tools listed in the table and note the according operating instructions of the dispersion tool.

Error codes

Any malfunctions during operation will be identified by an error message on the display (E and F).

Proceed as follows in such cases:

- ☛ Switch off device using the main switch (A)
- ☛ Carry out corrective measures
- ☛ Restart the device

Error code	Cause	Effect	Solution
Undervoltage			
Err	Mains voltage failed during operation.	Device off	- Switch device off and on
Err	Mains plug connected with switch set to "On".	Device off	- Switch device off and on
Overtemperature			
Err	Device will switch automatically to cooling mode if operated with excess load	Motor speed set automatically (not adjustable)	- Allow device to cool down: - Allow to continue running in cooling mode - Switch device off. Switch on again after some time. (Duration depends on the magnitude of the overtemperature) A standard power meter is used between the mains socket and the mains plug to prevent the device failing as a result of excess load. Power consumption < 500 W : continuous operation Power consumption > 500 W : device may switch to cooling mode
Electronic			
Err	Electronic fault	Motor speed set automatically (not adjustable)	- Switch off device Please send the device to the Servicing department for repair, enclosing a short description of the fault.

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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