37 846 00

IKA®

IKA ULTRA-TURRAX® T 50 basic

Betriebsanleitung	DE	5
Operating instructions	EN	15
Mode d'emploi	FR	25
Instrucciones de manejo	ES	35
Handleiding	NL	45
Istruzioni per l'uso	IT	55
Driftsanvisning	SV	65
Driftsinstruks	DA	74
Käyttöohje	NO	84
Driftsvejledning	FI	94
ΟΔΗΓΙΕΣ ΔΕΙΤΟΥΡΓΙΑΣ	EL	103
Instruções de serviço	PT	113

032009

Reg.-No. 4343-01

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¹ gilt für ≤ 100 mPa·s

Contents

	Page
Safety instructions	15
Proper use	16
Unpacking	16
Useful information	16
Drive unit	16
Assembly of boom and stand	17
Working with the unit	17
Dispersion tools: General notes	18
Dispersion tool S 50 N -	18
Dispersion tool S 50 KR -	19
Dispersion tool S 50 KG - HH -	19
Dispersion tool S 50 TT - W 40 TT	20
Agitator shaft R 50	21
Maintenance and cleaning	21
Accessories	22
Applicable standards and regulations	22
Guarantee	22
Technical data	23
Technical data: shafts	24
Technical data: rotors and stator	24
Spare parts list drive unit	124
Spare parts diagramm	125

Safety instructions

- The unit is not suitable for manual operation.
- Observe all pertinent safety instructions and directives, as well as work-safety and accident-prevention regulations for laboratory use.
- The unit must not be operated without a dispersion tool.
- Never run dispersion tools dry, as the gasket and bearings will be

destroyed if the tools are not cooled by the medium.

- The gasket and bearings are made of PTFE and rustproof steel; the following points should therefore be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

- Vibrations may loosen the hand screw on the flange of the drive unit. For your own safety check that the hand screw is secure and tighten it if necessary.
- Glass vessels must always be secured with a clamp to prevent them spinning. When working with flasks elastic couplers must be used to prevent fracture of the glass.
- Switch on the drive only at the lowest speed setting. Then gradually increase the speed to the desired level.
- When working with the R 50 agitator shaft the safety cage must not be removed. The agitator shaft must be installed in the **T 50 basic** drive unit such that the bent safety cage prevents any direct contact with the shaft from in front.
- Only dispersion and agitating tools approved by IKA may be used!
- Even in the case of repair work, the unit must only be opened by a qualified engineer. The mains plug must be pulled before opening the unit. The voltage-carrying components inside the unit may continue to carry a voltage for some time after the mains plug is pulled.
- Caution: Pull the mains plug before changing the dispersion heads!
- Caution: Continuous operation will cause the surface of the mounting flange for the dispersion tool to become hot!
- The ventilation slots on the drive must not be obstructed!
- The device does not start up again after a current failure. Switch device off and on.
- The **T 50 basic** may not be used to process inflammable or explosive media.

Proper use

The **ULTRA-TURRAX® T 50 basic** is a dispersion unit suitable for manufacturing emulsions and dispersions in combination with a dispersion tool. The drive unit must be operated on a stand.

Unpacking

The **ULTRA-TURRAX® T 50 basic** is supplied fully assembled. The dispersion tools are also supplied fully assembled. Unpack the unit and the dispersion tools carefully and check thoroughly for damage. It is important that any damage is recognized on unpacking, as an immediate damage report may have to be made (post, rail or shipping firm). As supplied, the package contains: one **ULTRA-TURRAX® T 50 basic** drive unit, one boom arm, one hexagon nut, one set of tools and one set of operating instructions.

Useful information

You have acquired a high-quality product. Used in combination with an S 50 series dispersion tool, the **ULTRA-TURRAX® T 50 basic** unit is a high-speed dispersing and emulsifying unit capable of handling free-flowing and liquid media in batches.

Dispersion is the dissolution and diffusion of a solid, liquid or gaseous phase in a continuum that is not consolute with that phase.


The rotor/stator principle:

Due to the high rotation speed of the rotor, the medium to be processed is automatically drawn axially into the dispersion head and then forced radially through the slots in the rotor/stator arrangement.

The high accelerations acting on the material produce extremely strong shear and thrust forces. In addition, high turbulence occurs in the shear gap between rotor and stator, which provides optimum mixing of the suspension.

The dispersion effectiveness is heavily dependent on the product of the shear gradient and the time the particles spend in the shear zone. The optimum range for the circumferential velocity of the rotor/stator arrangement is 10-24 m/s.

A processing time of a few minutes is usually sufficient to produce the desired fineness. Long processing

times bring only insignificant improvements in the obtainable fineness; the energy expended serves merely to increase the temperature of the medium.

Drive unit

With its power consumption. 1100 W at 10,000 rpm, the drive unit is suitable for a wide range of everyday dispersion applications in the laboratory.

The speed of the **T 50 basic** dispersion unit can be infinitely adjusted by using the numbered adjustment wheel.

An IKA-DZM control revolution counter can be inserted into the bore hole (10) in order to measure the speed. (cf. Fig. 1)

Assembly of boom and stand

Overview (Fig. 1)

The extension arm is mounted as follows:

- Insert the extension arm (2) into the flange
- Screw in cylinder bolt (1)
- Fasten cylinder bolt (1) using offset screwdriver SW4

Vibrations may loosen the hexagon nut. For your own safety check occasionally that the boom is securely fastened. Tighten the nut if necessary.

For safety the drive unit should be fastened to the telescopic stand R 2723 (4) or the stand R 2722 (4) using the cross sleeve R 261 (3). The weight compensation of the telescopic stand is of considerable benefit when raising and lowering the unit.

In order to increase the stability of the apparatus, the **T 50 basic** drive unit must be mounted as close as possible to the bar of the stand.

The dispersion tool (5) is pushed all the way into the drive unit and fastened with the hand screw (6). The distance between the dispersion tool and the base of the vessel should not be less than 20 mm.

The whole unit may be arranged slightly off-center to avoid an undesirable intake of air due to high rotation turbulence.

Working with the unit

Check that the voltage listed on the designation plate matches the mains voltage. Also pay attention to the ambient conditions listed in the technical data.

Only connect the device to the mains supply with the switch in position 0.

The driving unit is switched on by means of the On/Off- switch (pos. 6).

Error rectification

If the device doesn't run with the switch in position 1, switch the device off and back on again.

If these conditions are met the unit is ready for service when the mains plug has been plugged in. If this is not the case, safe operation cannot be guaranteed and the unit may be damaged.

Speed adjustment overview (Fig. 1)

Set the adjustment wheel (8) to the lowest speed before switching on the unit. The speed can be infinitely adjusted in accordance with the requirements of the medium to be processed by using the adjustment wheel. The no-load speed of the drive unit can be set to between 5,00 and 10,000 rpm.

Speed table

Note: The drive is fitted with control electronics which keep the speed constant for a wide range of applications. The speed can only be reduced when the upper limit has been reached (usually only with highly viscose media and large dispersing elements). You will need to connect a speed gauge (see "Accessories") if an exact speed setting and monitoring are required.

Adjusting wheel scale	1	2	3	4	5	6
Speed ref. to 30 l water	500	2000	4000	6000	8000	10.000

The speed table was calculated on the basis of a dispersion volume of 30 l water using the shaft S 50 N - G 45 M. The speeds depend on the dispersing element used as well as the viscosity and quantity of the medium. It may, therefore, change over the course of a dispersing process if the viscosity of the liquid changes.

Dispersion tools: General notes

Overview (Fig. 2), (Fig. 3), (Fig. 4)

Shaft:

The shaft mainly consists of the shaft tube, the rotary shaft, the bearing and the shaft seal. There are three types of shaft for the **ULTRA-**

TURRAX® T 50 basic:

S 50 N - / S 50 KR - / S 50 KG-HH -

and two additional agitating and mixing tools:

R 50 - Rührwelle / S 50 TT - W 40 TT

The shaft cannot be used alone. Only combinations of shaft and generator can be used for working with the medium. The shaft is supplied with a cap (remove before use). Use this cap to protect shaft tubes from damage when not in service.

Generator:

The generator (dispersion head) consists of a rotor and a stator. These parts can be unscrewed from the shaft. Note that rotors and stators form matching pairs and must not be confused. Nine different dispersion heads are available for various applications.:

- G 40 G / - G 45 G
- G 40 M / - G 45 M
- G 40 F / - G 45 F
- W 65 SK / - W 80 SMK

The dispersion tools you have ordered (e.g. S 50 N - G 45 G) are always supplied fully assembled. Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.

Dispersion tool S 50 N -

Overview (Fig. 4)

The bearing of the S 50 N - shaft is a friction bearing. This friction bearing

is made of PTFE. The shaft has a dispersion head added when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (21) and blade (22) from the rotary shaft (2001) using the generator and shaft wrenches included in the tool kit. The generator wrench is used to hold the rotor while unscrewing it by turning the shaft from above with the shaft wrench.

Caution: The dispersion head G 45 G has two prominent blades.

- Unscrew the stator (20) from the shaft (8) (left-hand thread). Insert the cross pin (3mm) of the shaft wrench into the slot on the end face of the shaft tube. The stator can now be unscrewed by turning it in the clockwise direction with the generator wrench.

- Pull the rotary shaft (2001) down and out of the shaft tube. The components (5, 10 and 11) threaded onto the shaft can simply be pulled off.

- The bearing bushing (5) must be replaced if it is worn.

Reassemble the dispersion tool in the reverse sequence. When assembling ensure that the bearing bushing and the rotor are correctly installed.

Spare parts :	Item	Designation
	2	Grooved pin
	5	PTFE bearing bushing
	8	Shaft tube
	10	Shaft sleeve
	11	Bearing bushing
	2001	Rotary shaft assembly

Sterilization:

After cleaning the dispersion tools the methods A), B) and C) are suitable for the S 50 N - dispersion tool. (See Maintenance, cleaning and sterilization of dispersion tools).

Risk of injury !

Risk of injury !

Dispersion tool S 50 KR -

Overview (Fig. 4)

The bearing of the S 50 KR - shaft consists of a ball bearing with FPM (FKM) radial shaft sealing rings. The shaft has a dispersion head added when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (21) and blade (22) from the rotary shaft (2001) using the generator and shaft wrenches included in the tool kit. The generator wrench is used to hold the rotor while unscrewing it by turning the shaft from above with the shaft wrench.

Caution: The dispersion head G 45 G has two prominent blades.

- Unscrew the stator (20) from the shaft (11) (left-hand thread). Insert the cross pin (3mm) of the shaft wrench into the slot on the end face of the shaft tube. The stator can now be unscrewed by turning it in the clockwise direction with the generator wrench.
- Pull the sealing bushing (13) out of the shaft (11). The radial shaft sealing rings can be pressed out of the sealing bushing using a suitable drift. When installing the radial shaft sealing rings the gap between the rings must be packed with grease.
- Press the rotary shaft (2001) with ball bearing (5) down and out of the shaft tube. The O-ring (7) will fall out of the shaft. Pull the shaft sleeve (9) off the rotary shaft (2001).
- The bearing bushing (9) must be replaced if it is worn.

Reassemble the dispersion tool in the reverse sequence.

Spare parts :	Item	Designation
	2	Grooved pin
	5	Deep-groove ball bushing
	7	O-ring
	8	Washer
	9	Shaft sleeve
	11	Shaft tube
	13	Sealing bushing
	2001	Rotary shaft assembly

Sterilization :

After cleaning the dispersion tools the method A) is suitable for this dispersion tool. The methods B) and C) must not be used otherwise the ball bearing and seals will be destroyed.

(See Maintenance, cleaning and sterilization of dispersion tools).

Dispersion tool S 50 KG - HH -

Overview (Fig. 4)

The bearing of the S 50 KG-HH - shaft consists of a ball bearing and a one-way friction ring seal. The friction rings are made of carbide. The shaft has a dispersion head added when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (21) and blade (22) from the rotary shaft (2001) using the generator and shaft wrenches included in the tool kit. The generator wrench is used to hold the rotor while unscrewing it by turning the shaft from above with the shaft wrench.

Caution: The dispersion head G 45 G has two prominent blades.

- Unscrew the stator (20) from the shaft (12) (left-hand thread). Insert the cross pin (3mm) of the shaft wrench into the slot on the end face of the shaft tube. The stator can now be unscrewed by

turning it in the clockwise direction with the generator wrench. Hold the dispersion tool with the stator vertical at the top to prevent the small springs of the friction ring seal from falling out.

- In this position the components of the friction ring seal can be pulled off the rotary shaft (2001).
- Press the rotary shaft (2001) with ball bearing (4) down and out of the shaft tube.
- The friction rings must be replaced if they are worn. The friction surfaces of the seal must be clean.

Reassemble the dispersion tool in the reverse sequence.

Spare parts:	Item	Designation
	2	Grooved pin
	4	Deep-groove ball bushing
	5	Circlip
	8	Friction ring seal assembly
	12	Shaft tube
	2001	Rotary shaft assembly
	3007	Spring
	3008	Flat seal
	3009	Round seal
	3010	Pin
	3011	Round seal

Sterilization :

After cleaning the dispersion tools the method A) is suitable for this dispersion tool. The methods B) and C) must not be used otherwise the ball bearing and the friction ring seal will be destroyed.
(See Maintenance, cleaning and sterilization of dispersion tools).

Dispersion tool S 50 TT - W 40 TT

Overview (Fig. 3)

This dispersion tool meets all requirements for corrosion-resistance. The parts coming into contact with the product are made from glass-reinforced PTFE and titanium.

The PTFE components (shaft tube, stator and bearing bushing) are relatively soft and thus permit operations in the liquid/liquid and liquid/gaseous ranges.

Because PTFE softens at high temperatures a maximum medium temperature of 100 °C must not be exceeded.

Hard, abradant and abrasive materials must not be processed.

High flow speeds (up to 20 m/s) result in chemical corrosion of the titanium components when used in high-oxidizing acids.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (15) from the rotary shaft (9) using the generator and shaft wrenches included in the tool kit. The generator wrench is used to hold the rotor while unscrewing it by turning the shaft from above with the shaft wrench.
- Unscrew the stator (14) from the shaft (10) (left-hand thread). Insert the cross pin (3mm) of the shaft wrench into the slot on the end face of the shaft tube. The stator can now be unscrewed by turning it in the clockwise direction with the generator wrench.
- Pull the rotary shaft (9) down and out of the shaft tube. The components (7; 11 und 12) threaded on to the shaft can simply be pulled off.
- The bearing bushing (7) must be replaced if it is worn.

Reassemble the dispersion tool in the reverse sequence. When assembling ensure that the bearing bushing is correctly installed. The stator must only be screwed on to the shaft tube by hand.

Risk of injury !

Risk of injury !

Spare parts :	Item	Designation
	2	Grooved pin
	7	PTFE bearing bushing
	8	Adapter
	10	Shaft tube
	12	Bushing
	13	Shaft sleeve
	14	Stator
	15	Rotor
	2001	Rotary shaft assembly

Agitator shaft R 50

Overview (Fig. 2)

The **R 50** agitator shaft is a high-speed ball-bearing rotary shaft. A splash guard seals it against the medium. The maximum immersion depth of the unit is 180 mm.

The free-wheeling conical shaft is dangerous. Observe the safety instructions !

The dissolver disc and cage can be removed for cleaning purposes.

Spare parts :	Item	Designation
	3	Grooved pin
	4	Deep-groove ball bearing
	5	Circlip
	6	Circlip
	7	Shim
	8	O-ring
	9	Bearing bushing
	10	Shim
	11	Rotary shaft
	13	Dissolver disc

14	Cage
16	Static ring
17	Friction ring
18	O-ring

Maintenance and cleaning

Drive unit:

The drive unit is maintenance-free but not immune from wear. The carbon brushes of the motor wear down over time and produce grooves in the collector. To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Dispersion tools:

The seals in the dispersion tools must be constantly monitored. In the event of leakage, the suction effect of the rotating shaft can cause the medium to penetrate as far as the drive unit. If liquid emerges from the side holes at the top of the shaft tube, stop work immediately and check the seals.

The function of the dispersion heads depends on the condition of the sharp edges on the rotor and stator. These edges may be blunted very quickly in abrasive media, reducing the effectiveness of dispersion.

The dispersion tools are cleaned by running them in a solvent which will dissolve substance residues without harming the seals. Because of the high flow speed, this is usually sufficient to clean the rotor and stator. The dispersion tool must be removed and cleaned immediately after use to prevent residues from adhering to the threads of the rotor and stator and to prevent the formation of undesirable bacteria cultures.

The dispersion tools may also be sterilized. The following methods may be used:

A) Chemical methods

Germicidal solutions (formalin, phenol, alcohol, etc.) are good general-purpose disinfectants. It is important that all disinfectant residues are then removed with sterilized water.

B) Sterilization with moist heat. This denotes the use of a steam jet pressurized to 2 bar at 120 °C.

C) Sterilization with hot air. Hot-air sterilization is usually performed at 160 °C or 190 °C.

When ordering spare parts, please state the fabrication number on the designation plate, the type of unit and the designation of the spare part.

Warranty

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Accessories

- DZM control Rpm counter
- Measuring cable
- R 2723 Telescopic stand
- R 271 Cross sleeve
- DK 50 Flow chamber
- R 2722 Stand
- R 261 Cross sleeve
- SV NS 60 Screw clamp

For R 50 agitator shaft:

- R 1402 Dissolver disk
- R 1405 Propeller
- R 1413 Turbine
- R 1414 Turbine

Applicable standards and regulations

EU Directives

EMC Directive: 89/336/EEC

Machines Directive: 89/392/EEC

EU standards

EN 61010-1:1993, VDE 0411-114

EN 50081:92, EN 55014:93, EN 50082:92, EN 60555:87

EN 292-1, EN 292-2, EN 60 204-1, EN 414

International standards

CAN/CSA C22.2 (1010.1:92), UL 3101-1

Technical data

Drive unit:

Speed range	rpm	500 - 10000
Speed display		Scale
Speed variation on load change	%	<1
Perm. ambient temperature	°C	5 - 40
Perm. humidity	%	80
Perm. on time (drive unit)	%	100
Protection to DIN EN 60529		IP20
Power consumption	W	1100
Power output	W	700
Design voltage	V AC	230 ±10%
	or	V AC
		115 ±10%
Frequency	Hz	50/60
Noise level	dBA	72 (without disp. tool)
Viscosity		See table, p.24
Dimensions (drive)	WxDxH	mm
Dimensions (boom)		mm
Weight		kg
Contamination level		2
Protection class		2
Overvoltage category	II	
Operation at a terrestrial altitude	m	max. 2000 above sea level

Subject to technical changes!

Technical data: shafts

Seal type	Shaft tube dia. mm	immersion depth mm	Materials in contact with product	pH range	Solvent	Abrasive substances ¹	Use up to °C ²	Vacuum to mbar	Pressure to MPa	Sterili- zation
S 50 N	27	70-250	1.4301 PTFE	2 - 13	yes	yes	180	-	-	yes, all methods
S 50 KR -	27	70-260	1.4301 FKM	2 - 13	no	no	80	100	-	yes, wet chemical
S 50 KG - HH -	27	70-260	1.4301 Carbide FKM	2 - 13	Conditional	no	120	30	-	yes, wet chemical
S 50 TT - W 40 TT	27	70-250	Titanium PTFE	1 - 14	yes	no	80	-	-	yes, wet chemical

¹ Wear should be expected above Mohs' hardness 6 () = Conditional
² Applies for medium temperature without additional heating by shear energy

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IKA® Korea Co LTD
IKA® Works Inc.