IKA ULTRA-TURRAX® T 18 basic

<i>BETRIEBSANLEITUNG</i>	<i>DE</i>	<i>4</i>
<i>OPERATING INSTRUCTIONS</i>	<i>EN</i>	<i>13</i>
<i>MODE D'EMPLOI</i>	<i>FR</i>	<i>22</i>
<i>INSTRUCCIONES DE MANEJO</i>	<i>ES</i>	<i>31</i>
<i>HANDLEIDING</i>	<i>NL</i>	<i>40</i>
<i>ISTRUZIONI PER L'USO</i>	<i>IT</i>	<i>49</i>
<i>DRIFTSANVISNING</i>	<i>SV</i>	<i>58</i>
<i>DRIFTSINSTRUKS</i>	<i>DA</i>	<i>67</i>
<i>DRIFTSVEJLEDNING</i>	<i>NO</i>	<i>76</i>
<i>KÄYTTÖOHJE</i>	<i>SF</i>	<i>85</i>
<i>INSTRUÇÕES DE SERVIÇO</i>	<i>PT</i>	<i>94</i>

Reg.-No. 4343-01

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Contents

EN

	Page
CE-Declaration of conformity	2
Safety instructions	13
Correct use	15
Unpacking	15
Useful information	16
Drive	16
Assembly of boom and stand	17
Working with the unit	17
Maintenance and cleaning	18
Accessories	18
Speed table	18
Permitted dispersion tools	19
Faults and troubleshooting	20
Warranty	21
Technical data	21

Safety instructions

For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - splashing liquids
 - body parts, hair, clothing and jewellery getting caught.


Fig. 3: Points of danger

- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
 - Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids etc.
 - **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
 - There may be electrostatic discharges between the medium and the dispersing instrument shaft which could pose a direct danger.
 - The appliance is not suitable for manual operation.
 - Set up the stand in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that the stand does not start to move.
 - The agitated vessels used for stirring have to be secured. Consider on a good stability of the entire structure.
 - Secure the stirring vessel against twisting.
 - Glass vessels must always be secured with a clamp to prevent them spinning. When working with flasks elastic couplers must be used to prevent fracture of the glass.
 - Note the operating instructions of the dispersing tool and accessories.
 - Only dispersing elements approved by **IKA** may be used.
 - Assemble exactly the dispersion tool following the instructions.
 - Please observe the permitted speed for the dispersing element. Never set higher speeds.
 - Do not use the appliance without a dispersing element.
 - Use the dispersing tool always inside the stirring vessel.
 - Check that the turning handles are secure and tighten if necessary.
 - Make certain that the unit is set at the lowest speed before commissioning; otherwise, the unit will begin running at the speed last set. Gradually increase the speed.

- Reduce the speed if the medium splashes out of the vessel because the speed is too high.
- Before switching on the dispersing instrument make sure that its shaft is immersed min. 45 mm in the medium to prevent the medium from splashing out.
- The distance between the dispersion tool and the vessel bottom should not be less than 10 mm (fig. 1).
- In the event of unbalance or unusual noises, switch off the appliance immediately. Replace the dispersing element. If there is no difference after the change of the dispersing tool, return it to the dealer or the manufacturer along with a description of the fault.
- Do not touch rotating parts during operation.
- Please note that the dispersing element and the journal bearings can become extremely hot during use.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do not process any flammable or combustible materials.
- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA** if you have any questions.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- Meter powder not too close to the flange. Powder can be blown away by air turbulences of the drive.
- The appliance starts up again following a cut in the power supply. Please note that the appliance must be switched off at the on/off switch in this case before the power supply is reconnected.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Always disconnect the plug before fitting accessories.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.
- The gasket and bearings are made of PTFE and rustproof steel;

the following points should therefore be noted: *Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.*

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

For protection of the equipment

- The voltage stated on the nameplate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts.
- Do not cover the ventilation slots on the drive in order to ensure adequate cooling of the drive.

Correct use

- **Use**

The **ULTRA-TURRAX® T 18 basic** are dispersion units which, in connection with a dispersion tool, manufacture emulsions and dispersions.

- **Range of use**

- Laboratories
- Pharmacies
- Schools

The disperser is conceived for the enterprise at a stand.

This device is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

Unpacking

- **Unpack**

- Please unpack the device carefully
- In the case of any damage a fact report must be set immediately (post, rail or forwarder)

- **Delivey scope**

- Drive unit **ULTRA-TURRAX® T 18 basic**
- Extension arm
- Hexagon socket screw
- Tool set
- Operating instructions

Useful information

Used in combination with an S 18-series dispersion tool, the drive unit is a high-speed dispersing and emulsifying unit capable of handling free-flowing and liquid media in batches.

Dispersion is the dissolution and diffusion of a solid, liquid or gaseous phase in a liquid that is not consolute with that phase.

The rotor/stator principle

Due to the high rotation speed of the rotor, the medium to be processed is automatically drawn axially into the dispersion head and then forced radially through the slots in the rotor/stator arrangement. The high accelerations acting on the material produce extremely strong shear and thrust forces. In addition, high turbulence occurs in the shear gap between rotor and stator, which provides optimum mixing of the suspension.

The dispersion effectiveness is heavily dependent on the product of the shear gradient and the time the particles spend in the shear zone. The optimum range for the circumferential velocity of the rotor/stator arrangement is 6-24 m/s.

A processing time of a few minutes is usually sufficient to produce the desired fineness. Long processing times bring only insignificant improvements in the obtainable fineness; the energy expended serves merely to increase the temperature of the medium.


Fig. 4

Drive unit

The drive unit opens up a wide range of possibilities for dispersion technology in the conventional laboratory with a respective performance output of approximately 300 watts at 24,000 rpm.

The speed of the **ULTRA-TURRAX® T 18 basic** dispersion unit is set by using the differently adjusting wheel. In order to measure the speed of the dispersion unit an **IKA DZM** control tachometer may be connected after removal of the cover (pos. 3), (see fig. 1).

Assembly of boom and stand

The delivered extension arm is mounted as follows (overview fig. 1):

- Insert extension arm (pos. 1) into flange
- Screw in the cylinder head screw (pos. 2)
- Firmly tighten cylinder head screw (pos. 2) using the A/f 4 right angled screwdriver.

Vibrations may loosen the cylinder head screw. Therefore, for your own safety, periodically check that the extension arm is secure. If necessary, retighten the cylinder head screw.

To ensure safe operation the driving unit is fastened to the plate stand R 1826 (pos. 8) by means of a cross sleeve (pos. 7).

In order to enhance the stability of the mechanical arrangement, the driving unit has to be mounted as close as possible to the stand pillar.

Working with the unit

Assembling of the dispersion tool (overview fig. 1)

- Open the knurled screw (pos. 10) so that the thread does not protrude into the hole.
- Insert the dispersing element (pos. 9) from the bottom into the drive unit as far as the stop. Ensure that the dispersion tool is aligned correctly when inserted (refer to picture: the tabs must project sideways when viewed from the front of the drive unit/dispersion tool). The shaft will audibly click into place when a small amount of force is applied. The dispersion tool is correctly mounted when the edge of the projecting part is aligned with the lower edge of the drive unit (see fig. 5).
- Now secure the dispersing element by screwing the knurled screw tight (pos. 10).



Fig. 5

The filling level of the medium may not be less than about 55 mm.

The entire unit may also be arranged somewhat off centre in order to avoid any undesired air admission caused by the strong rotational turbulence.

The adjustment wheel (pos. 4) is set at the lowest speed before the unit is switched on. The speed which is set may be read by comparing the colors on the scale (pos. 4) and the adjustment wheel.

By using the adjustment wheel (pos. 4) the speed may be infinitely varied to meet the requirements of the medium to be processed. The no-load speed of the **ULTRA-TURRAX® T 18 basic** driving unit may be set between 3,500 and 24,000 rpm.

If these conditions are met the unit is ready for service when the mains plug has been plugged in.

The driving unit is switched on by means of the On/Off- switch (pos. 6).

Maintenance and cleaning

Drive unit: The drive unit is maintenance-free but not immune from wear. The carbon brushes of the motor wear down over time.

For cleaning disconnect the main plug.

To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Dispersion tools: Consider please the manuals associated to the dispersion tools.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see **www.ika.de**, spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Accessories

- DZM control.o Rpm counter
- R 1826 Plate stand
- R 182 Cross sleeve
- RH 3 Tensioner

Speed table

Note: the drive is fitted with control electronics which keep the speed constant for a wide range of applications. The speed can only be reduced when the upper limit has been reached (usually only with highly viscose media and large dispersing elements). You will need to connect a speed gauge (see "Accessories") if an exact speed setting and monitoring are required.

Adjusting wheel scale	1	2	3	4	5	6
No-load speed	3.600	7.200	11.200	15.600	20.000	24.200
Speed with 0,1l water (rpm)	3.500	7.000	11.000	15.500	20.000	24.000

The speed table was calculated on the basis of a dispersion volume of 0.1l water using the dispersing element S18N-10G (at nominal tension 230V respectively 115V). The speeds depend on the dispersing element used as well as the viscosity and quantity of the medium. It may, therefore, change over the course of a dispersing process if the viscosity of the liquid changes.

Permitted dispersion tools

Designation		Schaftwerkstoff
S18N-10G		Stainless steel
S18N-19G		Stainless steel
S18D-10G-KS		Synthetic material
S18D-14G-KS		Synthetic material

For applications and further information, please see the operating instructions of the dispersing element.

Only use the dispersion tools listed in the table and note the according operating instructions of the dispersion tool.

Faults and Troubleshooting

Please follow the following procedure should a fault occur during operation:

- ☞ Switch off device using the main switch (6)
- ☞ Carry out corrective measures
- ☞ Restart the device

Fault	Effect	Solution
Undervoltage		
Mains voltage failed during operation	Device off	- Switch device off and on
Mains plug connected with switch set to "On"	Device off	- Switch device off and on
Overtemperature		
Device will switch automatically to cooling mode if operated with excess load	Motor speed set automatically (not adjustable)	- Allow device to cool down: - Allow to continue running in cooling mode - Switch device off. Switch on again after some time. <i>(Duration depends on the magnitude of the overtemperature)</i> A standard power meter is used between the mains socket and the mains plug to prevent the device failing as a result of excess load. Power consumption < 500 W : continuous operation Power consumption > 500 W : device may switch to cooling mode
Electronic fault		
Device switches to cooling mode	Motor speed set automatically (not adjustable)	- Switch off device and allow to cool in order to exclude the possibility of overtemperature. If the device continues to switch to the cooling mode after it has cooled down then please send it to the Servicing department for repair, enclosing a short description of the fault.

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Technical data

Speed range	rpm	3,500 - 24,000 (at nominal-voltage 230 V and 50 Hz)
Speed display		Scale
Speed variation on load change	%	<6
Perm. ambient temperature	°C	5 - 40
Perm. humidity	%	80
Perm. On- time (drive unit)	%	100
Protection to DIN EN 60529		IP20
Power consumption	W	500
Power output	W	300
Design voltage	VAC	230 ±10%
or	VAC	115 ±10%
Nominal voltage	VAC	230
or	VAC	115
Frequency	Hz	50/60
Noise level (without disp. tool)	dB(A)	73
Dimensions (drive) W x D x H	mm	65 x 80 x 240
Dimensions (boom)	mm	Ø13 / L160
Weight	kg	1,6
Contamination level		2
Protection class		II
Overvoltage category		II
Operation at a terrestrial altitude	m	max. 2000 above sea level

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