IKA® LR-2.ST

<i>BETRIEBSANLEITUNG</i>	<i>D</i>	<i>3</i>
<i>OPERATING INSTRUCTIONS</i>	<i>GB</i>	<i>14</i>
<i>MODE D'EMPLOI</i>	<i>F</i>	<i>25</i>
<i>INDICACIONES DE SEGURIDAD</i>	<i>E</i>	<i>37</i>
<i>VEILIGHEIDSINSTRUCTIES</i>	<i>NL</i>	<i>38</i>
<i>NORMI DI SICUREZZA</i>	<i>I</i>	<i>39</i>
<i>SÄKERHETSANVISNINGAR</i>	<i>S</i>	<i>40</i>
<i>SIKKERHEDSHENVISNINGER</i>	<i>DK</i>	<i>41</i>
<i>SIKKERHEDSHENVISNINGER</i>	<i>N</i>	<i>42</i>
<i>TURVALLISUUSOHJEET</i>	<i>SF</i>	<i>43</i>
<i>Υποδείξεις ασφάλειας</i>	<i>GR</i>	<i>44</i>
<i>INSTRUÇÕES DE SEGURANÇA</i>	<i>P</i>	<i>45</i>

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D

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CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60010-1; EN 50082; EN 55014; EN 55114 and EN 60555.

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F

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60010-1; EN 50082; EN 55014; EN 55114 and EN 60555.

IKA WERKE GmbH & CO.KG

Staufen, Mai 2003

Reiner Dietsche
Vice President

Wolfgang Buchmann
Corporate Director - Quality Assurance

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Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Explanation of symbols

This symbol identifies information that is absolutely essential for your safety and health. Failure to observe this information may result in negative effects for your health and injuries.

This symbol indicates information which is necessary to ensure the safe handling and operation of the laboratory reactor. Failure to observe this information could result in incorrect test series or damage to the appliance.

Contents

	Page
Guarantee	2
Explanation of symbols	14
Contents	14
Safety instructions	14
Proper use	15
Unpacking	16
System configurations	16
Setting up and operation	16
LR-2.ST stand system	17
Stirrer	18
Reactor head	18
Reaktor vessel and clamping mechanism	19
Safety disconnection LR-2.SI	19
Splinter guard LR-2.SP	20
Use of dispersing element	21
Cleaning / sterilization	21
IKA® Visco-modul VM 600 basic	22
Maintenance	22
Accessories	23
Materials in contact with medium / laboratory reactor	24
Technical data	24
List of spare parts and spare parts diagram	47 - 51

Safety instructions

You must use suitable personal protective equipment (PPE) when using the laboratory reactor system!
All work with the reactor system must be carried out behind a protective screen or under a closed fume hood!

An EMERGENCY STOP switch must be fitted between the laboratory reactor and the mains power supply. This must completely disconnect the equipment from the power supply in the event of an emergency. You must use the reactor with the safety disconnector and the splinter guard provided.

Connect the laboratory stirrer and the disperser drive to the mains power supply over the LR-2.SI safety disconnector. Connect these appliances to the multiset on the safety disconnector only.

Appliances that need to continue operation once the reactor cover has been opened (e.g. thermostat), must not be operated via the safety disconnector.

Implement the accident prevention regulations for this type of activity and workplace. Wear your personal protective equipment at all times when working with the laboratory reactor.

Do not use flammable or explosive gases, vapours or solvents in the LR-2.ST laboratory reactor!
The appliance is not designed for use in explosion hazard (Ex-certificate) areas!

The equipment is not designed for overpressure use.

The appliance is designed for operation at vacuum settings up to 25 mbar (see Accessories).

Please follow the safe procedures specified for handling reactants and reaction products.

Certain applications and materials may be hazardous. You should take precautions to prevent contact with, or inhalation of, toxic liquids, gases, fumes, vapours or powders.

Risks may also be posed by biological or microbiological substances.

Moving and rotating equipment parts also constitute a hazard.

The laboratory stirrer/disperser is fitted with a limit switch which prevents the unit from being started until it is in the fully lowered position. The anchor stirrer and the dispersing element must only be operated when the reactor vessel is fully closed.

Do not start up the laboratory reactor when it is open.

If the anchor stirrer or dispersing element rotate when they are covered in a product and the reactor is open, this will cause parts or liquids to be projected outwards.

The LR 2000.1 and LR 2000.2 laboratory reactor vessels may only be heated without pressure using the double shell with thermostats or similar heat sources. Do not under any circumstances use a heating mantle or a hot plate (Bursting hazard!).

We recommend that you use the IKA HBR 4 digital heating bath to temper the LR 2.1 single-walled laboratory reactor vessel.

The reactor system must always be ventilated when working under normal pressure in order to prevent any pressure build-up caused by highly volatile gases or unpredictable reaction pressure gradients. Condense volatile gases using a cooler with a standard ground connector (e.g. a reflux condenser) on the reactor cover.

The maximum permitted temperature inside the reactor vessel is 230°C. Do not exceed this temperature.

Burns Hazard. The heated reactor vessel and cover are very hot. Wear protective gloves when handling hot parts.

Use the splinter guard provided. It will also prevent you from accidentally touching the hot reactor vessel.

Ensure that the thermostat used for tempering is fully efficient. A defective thermostat can cause uncontrolled reactions.

Before you fill the reactor vessel, ensure that the reagents used do not corrode the FFPM seal.

Only use IKA® approved accessories!

Use only original IKA® spare parts!

Proper use

The IKA® LR-2.ST system is a modular miniplant reactor system. It has been designed to simulate and optimise chemical reaction processes and can also be used for mixing, dispersion and homogenisation in modelling processes.

The working volume is 500 to 2000 ml.

Thanks to the highly resistant FFPM seals the medium in the reactor vessel can be heated to 230°C. Vacuum operation is possible up to 25 mbar (see Accessories).

By adapting the VM 600 Basic Visco Module (see Accessories) you will also be able to use torque measurements to test the rheological properties and the flow and deformation behaviour of materials.

The IKA® laboratory equipment range is vast. It provides numerous opportunities for building equipment combinations customised for a wide variety of individual applications.

You can use these in conjunction with labworldsoft® laboratory software to create a high-performance system for measuring, controlling, regulating and fully documenting all laboratory processes.

Sample applications:

- Manufacture of cremes, lotions, emulsions and liposome preparations in the pharmaceutical and cosmetic sector (lipstick, sun cream, eyeliner).
- Mixing of solids such as calcium carbonate, talc and titanium oxide with liquid polymers.
- Mixing of additives and solid polymer compounds with mineral oils.
- Dispersion and mixing of solids and fibres in liquids and polymers.
- Manufacture of milk products and dietary products.
- Extraction of enzymes from biomass.

Unpacking

The LR-2.ST laboratory reactor system is supplied as a fully functional unit in various configurations depending on customer requirements. Accessories, such as ULTRA-TURRAX, temperature gauges and reactor vessels are supplied separately.

Please unpack the individual components and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

System configurations

The basic IKA® Miniplant laboratory reactor comprises:

- LR-2.ST stand system complete with reactor cover
- EUROSTAR power control visc P7 drive unit (or EUROSTAR power control visc P4)
- LR 2000.11 anchor stirrer with flow borings
- LR-2.SP splinter guard
- LR-2.SI safety disconnecter

This system can be used with a **variety of reactor vessels**:

- LR 2.1 single-walled reactor vessel
- LR 2000.1 double-walled reactor vessel
- LR 2000.2 double-walled reactor vessel with bottom discharge valve

All reactor vessels are made of borosilicate glass 3.3 and have a volume of 2000 ml. The nominal width is DN150. The double-walled reactor vessels are all fitted with stainless steel quick-release connectors for thermostat hoses.

There are two tempering units for **tempering** the reactor vessels:

- HBR 4 digital heating bath (for single-walled reactor vessels)
- LT 6 control circulation thermostat (for double-walled reactor vessels)

Connect the thermostats and double-walled reactor vessels using two LT 5.20 hoses and LT 5.24 locking unions.

Another system configuratin for **testing rheological material properties** is the basic system plus

- IKA® - visco module VM 600 basic

Note: The VM 600 Basic Visco Module can only be used in vacuum processes and with the LR 2000.10 anchor stirrer with a PTFE scraper in a very limited set of conditions.

Other possible system components are listed in Accessories.

Setting up and operation

Follow the operating instructions given in this manual!

You should also refer to the operating instruction manual provided with the other items of equipment delivered to you!

The reactor system must be operated under a closed fume hood or behind a protective screen!

The reactor must only be used with hte safety disconnecter supplied. The safety disconnecter ensures that the system can be completely disconnectet from the mains power supply when required!

Always use the IKA® LR-2.SP splinter guard!

Only connect the laboratory stirrer and the disperser drive to the multi-socket. The multi-socket is mounted on the back of the equipment and is connected to the safety disconnecter.

Fit the safety disconnecter and the splinter guard!

Secure the reactor vessel using the lock ring on the reactor cover.

Connect the power supply cables to the mains supply.

Switch on the safety disconnecter by pressing the green button.

Start the laboratory stirrer.

The IKA®LR-2.SI has an EMERGENCY STOP button and also a safety disconnecter which provides an alternative means of shutting down the power supply circuit when the reactor cover is open. This protects you from rotating elements such as anchor stirrers or dispersing elements.

**The appliance will not restart automatically after a power failure.
You must switch the safety disconnecter back on.**

On system configurations which feature a single-walled reactor vessel, the HBR 4 digital heating bath also serves as a splinter guard. The glass reactor vessel must always be fully immersed in the tempering unit. With this equipment configuration you must ensure that any vapours released from the tempering bath are extracted.

Any reaction vapour released at the reactor cover must be extracted!

LR-2.ST stand system


The stand system holds the drive unit complete with the reactor cover and the stirring instrument, the reactor vessel and the safety disconnecter.

An adjustment screw on the rear, right-hand side foot of the stand can be used to level the stand and compensate for any unevenness in the mounting surface.

Undo the lock nut. Turn the knurled screw clockwise to slightly raise the foot of the appliance, turn it anticlockwise to lower it. When the appliance is level, tighten the lock nut.

If you are working with the LR 2.1 single-walled reactor vessel and the HBR 4 digital heating bath, you must remove the baseplate including clamping bracket and one side support from the stand.

Undo the top two M8 cheese-head screws using the DIN 911-6 cranked wrench key provided. Remove the baseplate, undo the bottom cheese-head screws and remove the side support.

Adjuster screw for levelling

M8 cheese-head screw with baseplate and side support

The remaining side support is used to fit the safety disconnecter. Place the plastic cover onto the aluminium profile.

To raise the whole assembly by hand, first undo the two knurled screws on the left and right of the stand legs.

The pressure springs built into the telescopic stand mean that the whole reactor can be easily raised and lowered.

Undo the knurled screw to adjust the height

The laboratory stirrer mount with integrated stirrer shaft bearing is fitted to the stand adapter (connected to the stand legs) with two tie rods. The ball-bearing, vacuum-tight stirrer shaft is connected to the stirrer shaft drive (EUROSTAR power control visc P7) by a

coupling flange. The stirrer shaft features covered precision bearings and is sealed with rotary shaft seals with PTFE sealing lips for vacuum operation.

Penetration of the stirrer shaft inevitably produces wear debris from the stirrer shaft seal. In order to keep this debris away from the product, a protective ring has been fitted to the stirrer shaft from below.

The anchor stirrer is also inserted into the stirrer shaft from below. The torque is transferred by a square shaft to the anchor stirrer. This is secured axially by means of a retaining ring.

Laboratory stirrer mount with stirrer shaft bearing

Temperature sensor (accessory)

Lock ring

Protective ring

Flow breaker (accessory)

Anchor stirrer

Stirrer

Refer to the operating instructions for the laboratory stirrer.


The laboratory stirrer is a fixed component of the laboratory reactor system.
The following IKA® laboratory stirrers are available for this system:

- EUROSTAR power control visc P7
 - EUROSTAR power control visc P4
- (or)

The EUROSTAR power control visc P4 is ideal for a maximum viscosity of 100 000 mPas.

The standard EUROSTAR power control visc P7 is used for a maximum viscosity of 150 000 mPas. (see Technical Data)

Reaktor head

The reactor cover is screwed onto the stirrer shaft bearing using a thread. The reactor cover can be easily removed for cleaning or sterilisation purposes (see Cleaning / Sterilisation).

An O-ring provides the seal between the stirrer shaft bearing and the reactor cover.

The reactor cover is fitted with three NS 29/32 standard ground connectors and two NS 14/23 standard ground connectors. These connectors are screwed to the cover and sealed by O-rings.

One of the NS 14/23 connectors is fitted with a ventilation and air vent valve, whilst the other connectors are closed off with DIN 12252 ball valves.

CAUTION:

Hazardous pressure build-up!

Ensure that the reactor vessel is properly vented and can release sufficient air!

Adjuster screw

All sealing rings on the reactor cover are made of FFPM (perfluoroelastomer).

Dispersing elements (ULTRA TURRAX), temperature sensors, flow breakers and other accessories can be fitted to the free connectors. In order to do this, unscrew the threaded connector plug from the cover using the pin provided and then fit the device (e.g. LR 2000.60 temperature sensor) in its place.

The lock ring used to secure the reactor is located on the rim of the cover.

The clamping force of the lock ring can be set when it is open using the adjuster screw.

Reactor vessel and clamping mechanism

Configuration with single-walled reactor vessel

The single-walled reactor vessel is directly centred with and connected to the reactor cover using the lock ring.

Position the stainless steel cover plate on the tempering bath and centre the vessel with the cork ring underneath the reactor cover. Lower the stand as far as the top glass edge and secure the glass vessel by closing the clamping lever on the lock ring.

When you insert the vessel, check that the sealing contact surface of the glass vessel is undamaged and that the seal on the reactor cover fits properly.


Configuration with double-walled reactor vessel

The LR 2000.1 double-walled reactor vessels without bottom discharge valve and LR 2000.2 double-walled reactor vessels with bottom discharge valve are set onto the cork ring provided.

Secure the reactor vessel using the clamping strap on the stand. Screw the LT 5.24 hose adapters to the LT 5.20 tempering hoses.

Quick-release connector

Then connected the fully assembled hoses to the LT 6 control circulation thermostat.

Lower the drive unit including reactor cover and anchor stirrer and connect the reactor vessel to the cover by closing the lock ring. Firstly check that the sealing surface on the glass vessel is intact and that the seal on the reactor cover fits properly.


Clamping strap

Double-walled vessel

Safety disconnection LR-2.SI

You can mount the safety disconnector power box on either the left or the right aluminium profile using two slider blocks. Insert the slider blocks into the profile groove with the spring-loaded ball to the back. You do not need to remove the baseplate. Position the slider blocks in the profile in such a way that you can screw on the power box with the two M6 cheese-head screws using the SW5 Allen key.

Slide the wireless safety switch down to the bottom end of the top (thicker) stand rod and secure it using the clamping screw.

The assembled safety switch should be pointing downwards. Turn the bracket so that it is positioned in such a way that it will not prevent you from operating the clamping screw for the stand properly.

Lower the stand until the reactor cover is sitting on top of the reactor vessel and secure it in this position by tightening the knurled screws.

Fit the second safety switch (with cable) to the bottom (thinner) stand rod in this position and secure it so that it is approx. 5 mm away from the top safety switch. Align the two switches with each other.

Attach the socket bracket to the top of the stand and place the multi-socket on top

Connect the mains plugs of the stirrer drive and the disperser drive (if you are using one) to the multi-socket on the safety disconnecter.

Fix the power cables together with cable ties. When working with the IKA®Visco module VM 600 basic it is particularly important that the power cable of the stirrer drive is secured to prevent a restoring torque.

Connect the power supply cable of the safety disconnecter to your power supply.

Switch on the laboratory stirrer drive at its main switch. Close the reactor cover and press the green button on the safety box. The stirrer drive will start up.

To switch off the stirrer drive, press the red button or open the reactor cover.

The IKA®LR-2.SI safety disconnecter disconnects the laboratory reactor system from the mains power supply when the red EMERGENCY STOP button is pressed.

The safety disconnecter must be actively switched on again if the power supply is interrupted. The appliances cannot automatically restart.

The limit switch trigger on the stand rod (switch distance approx. 5 mm) also breaks the circuit if the reactor cover is open. This prevents any danger from rotating elements such as anchor stirrers or dispersing elements.

The safety disconnecter must be switched back on again each time it is switched off (use the green button). After a power failure it will also be necessary to switch on the safety disconnecter.

Splinter guard LR-2.SP

The Macrolon guard is mounted on the baseplate and protects the operator from splinters and spray. It reduces the risk of injury caused by broken glass, and burns as a result of accidentally touching the hot reactor vessel.

The splinter guard can be quickly fitted and removed again without the need for special tools.

Place the Macrolon guard panel onto the baseplate and secure it using the two knurled screws on the bottom of the baseplate.

*Fixing the guard
in position*

CAUTION:

The guard does not provide protection against bursting!

DO NOT operate the LR-2.ST laboratory reactor at excessive pressures even when the guard is fitted.

Use of a dispersing element is only permitted when the reactor contains a working volume of at least 800 ml.

Using a dispersing element with volumes lower than this will result in dry running which will damage the axial face seal in the dispersing element.

The following dispersing elements can be used with the ULTRA TURRAX T 25 BASIC drive unit:

	S 25 KV - 18 G	S 25 KV - 25 G	S 25 KV - 25F
Rotor/stator gap width (mm)	0,3	0,5	0,5
Max. temperature (°C)	220	220	220
Suspension fineness (µm)	10-50	15-50	5-25
Emulsion fineness (µm)	1-10	1-10	1-5

All parts in contact with the medium are made of 1.4571 stainless steel or FFPM/SiC. The minimum/maximum immersion depth is 40/225 mm. These dispersing elements can be used in the reactor system at vacuums of up to 25 mbar.

The LR 2000.40 shaft holder is necessary to install the dispersing element in the reactor cover.

Circlip

Clamping screw

Circlip

Assembly:

Insert the circlip into groove on the shaft

Place shaft into the shaft holder from below

Insert the circlip into upper groove

Tighten clamping screw

The circlips required are supplied with the laboratory reactor.

The assembled shaft holder is replaced by a standard ground connector in the reactor cover.

Cleaning / sterilization

Cleaning

When the laboratory reactor is open the reactor vessel can be completely removed from the reactor so that the interior can be cleaned.

The axially inserted anchor stirrer can be removed from the stir shaft with a quick downward pull. The PTFE scrapers (LR 2000.10 anchor stirrer) can be removed from the anchor stirrer for thorough cleaning.

If no additional devices such as ULTRA TURRAX, temperature sensors or flow breakers are installed, the reactor cover can be unscrewed from the stirrer shaft bearing for cleaning.

First remove the connector components, the connectors for the additional devices and pull the protective ring off the stir shaft on the cover.

The stirrer shaft and the protective ring can also be cleaned now.

All O-rings can now be carefully removed by hand or using a blunt tool.

Removing the anchor stirrer

Chemical methods:

Germicidal solutions (formalin, phenol, alcohol etc.) are good general-purpose disinfectants. It is important that all disinfectant residues are then removed with sterilized water.

Sterilization with moist heat:

This denotes the use of a steam jet pressurized to 2 bar at 120 °C.

Sterilization with hot air:

Hot-air sterilization is usually performed at 160 °C or 190 °C.


Laboratory reactor components may only be sterilised using these methods.

IKA® Visco-modul VM 600 basic

Refer to the VK 600 operating instructions!

The IKA®Miniplant LR-2.ST laboratory reactor system can be fitted with a tool to measure the restoring torque. This can be used to measure rheological properties, flow and deformation behaviour.

The LR-2.ST laboratory reactor can also be combined with the VM 600 basic IKA® visco module at a later stage. In this case, the device should be shipped back to IKA® for assembly.

The module will be mounted between the laboratory stirrer drive and the stirrer shaft bearing.

Measurement range: 0 - 600 Ncm
 Reproducibility: +/- 0.1 Ncm static
 +/- 0,5 Ncm dynamic

RS-232 interface with analogue output
 offset compensation for error elimination

*Visco Modul
 VM 600 basic*

Note:

The VM 600 Basic Visco Module can only be used in vacuum processes and with the LR 2000.10 anchor stirrer with a PTFE scraper in a very limited set of conditions. This is because the friction losses are highly variable.

Maintenance

The IKA® laboratory reactor system requires little maintenance. The ball bearings are sealed and lubricated for life.

Refer to the operating instructions for individual components.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name and item-no. of the spare part.

Before you send any equipment to us for repair, ensure that it has been cleaned of any toxic or corrosive substances. Use only water with a detergent additive containing a surface active agent or, in cases of extreme contamination, isopropyl alcohol for cleaning parts not in contact with the medium.

Accessories

Range	Name	Designation	Remark
Mixing	LR 2000.10	Anchor stirrer with PTFE scraper	Limited torque measurement possible
	LR 2000.11	Anchor stirrer with flow borings	
	LR 2000.20	Flow breaker	
Tempering	HBR 4 digital	Heating bath	for LR2.1 single-walled reactor vessel
	LT 6 control	Circulation thermostat	for LR 2000.1 and LR 2000.2 double-walled reactor vessel
	LT 5.20	Hoses	for LT 6 control
	LT 5.24	Hose adapter M16x1	in connection with LT 5.20
Temp. measurement pH measurement	DTM 12	Digital temperature gauge or pH measuring instrument	with RS 232 interface
	PT 100.25	Temperature sensor	
	LR 2000.60	Temperature sensor holder	FFPM - solvent-based seal; maximum temperature 230 °C
	pH 225.1	pH elektrode with NTC temperature sensor	
	DTM 12.15	Cable for pH electrode	
	LR 2000.65	pH sensor holder	
Software	labworldsoft		Laboratory instrument software for control and data collection
Vacuum	IKAVAC VC2	Vacuumcontroller	to create a regulated vacuum
	MZ-2C	Chemistry - diaphragm pump	
	VC 2.4	Pump control	for use with an electric vacuum pump only in conjunction with VC2
Dispersing	ULTRA-TURRAX T 25 basic	Dispersing unit	
	S 25 KV - 18 G	Dispersing tool	
	S 25 KV - 25 G	Dispersing tool	
	S 25 KV - 25 F	Dispersing tool	
	LR 2000.40	Shaft holder	for Dispersing tool to 230 °C / FFPM seal
Reactor vessels	LR 2.1	Single-walled reactor vessel	Borosilicate glass 3.3 volumen 2000ml
	LR 2000.1	Double-walled reactor vessel	Borosilicate glass 3.3 volumen 2000ml
	LR 2000.2	Double-walled reactor vessel with bottom discharge valve	Borosilicate glass 3.3 volumen 2000ml
Safety	LR-2.SP	Splinter guard	Macrolon
	LR-2.SI	Safety disconnection	

Materials in contact with medium/laboratory reactor

Reactor cover	Stainless steel 1.4571
Reactor vessel	Borosilicate glass 3.3
Anchor stirrer	Stainless steel 1.4571
Scraper	PTFE
O - ring	FFPM
Shaft seal	PTFE / 1.4571

Technical data

Stirrer drive

Data	Unit	EUROSTAR power control visc P4	EUROSTAR power control visc P7
Speed range	1 / min	14 - 530	8 - 290
Maximum torque stirrer shaft	Ncm	200	380
Temporary overload	%	200	200
Speed display		digital	digital
Maximum power stirrer shaft	W	95	95
Overall efficiency		0,55	0,55
Viscosity range of application	mPas	to 100 000	to 150 000

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Data	Unit	
Design voltage	VAC	230 ±10%
	VAC	115 ±10%
Frequency	Hz	50 / 60
Perm. ambient temperature	°C	+5 to +40
Perm. relative humidity	%	80
Dimensions B x H x T (closed)	mm	460 x 430 x 900
Dimensions B x H x T (openend)	mm	460 x 430 x 1130
Weight (without reactor vessel)	kg	25
Test volume	ml	500 - 2000
Max. Temperature of tempering bath	°C	230
Pression de service	mbar	25 - 1013
Operation at a terrestrial altitude	m	max. 2000 above sea level

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

Please take spare parts of the optional system components from the instruction manuals!

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



IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.

2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd

No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.

293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.

Estrada do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-9600
FAX +55 21 2435-9601
E-mail: brasil@ika.net

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou

173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn