## RET control/t

Betriebsanleitung	DE	3
Operating instructions	EN	26
Mode d'emploi	FR	49
Indicaciones de seguridad	ES	72
Veiligheidsinstructies	NL	74
Norme di sicurezza	IT	76
Säkerhetsanvisningar	SV	78
Sikkerhedshenvisninger	DA	80
Sikkerhetshenvisninger	NO	82
Turvallisuusohjeet	FI	84
Instruções de segurança	PT	86
Wskazówki bezpieczeństwa	PL	88
Bezpečnostní upozornění	CS	90
Biztonsági utasítások	HU	92
Varnostna navodila	SL	94
Bezpečnostné pokyny	SK	96
Ohutusjuhised	ET	98
Drošības norādes	LV	100
Saugos reikalavimai	LT	102

022009

Reg.-No. 4343-01

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DE

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## EG-DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EC and 2006/95/EC and conforms with the standards or standardized documents DIN EN IEC 61010-1, -2-010 and DIN EN IEC 61326-1.

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In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

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## 1 Warranty

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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# Contents

	Page
<b>1 Warranty</b>	<b>2</b>
<b>2 Safety instructions</b>	<b>27</b>
<b>3 Correct use</b>	<b>28</b>
<b>4 Unpacking</b>	<b>28</b>
<b>5 Commissioning</b>	<b>29</b>
<b>6 Control elements</b>	<b>29</b>
<b>7 Display</b>	<b>29</b>
<b>8 Switching on</b>	<b>30</b>
<b>9 Factory settings</b>	<b>30</b>
9.1 Resetting parameters to factory settings (RESET)	31
<b>10 Operation modes</b>	<b>31</b>
10.1 Operation mode A	31
10.2 Operation mode B	31
10.3 Operation mode C	31
10.4 Operation mode X1	31
10.5 Changing the operation modes	31
<b>11 Function “Heating”</b>	<b>32</b>
11.1 Medium temperature control by temperature sensor PT 100	32
11.2 Medium temperature control by twin temperature sensors PT 1000	32
11.3 Medium temperature control without external temperature sensors	33
11.4 Residual heat warning (HOT)	33
11.5 Changing the temperature unit of measurement from °C; °F	34
11.6 Setting the safety temperature	34
<b>12 Function “Stirring”</b>	<b>34</b>
12.1 Intermittent mode	34
12.2 Stirring bar breakaway monitoring	35
12.3 Setting the intermittent mode and the stirring bar breakaway monitoring function	35

	Page
<b>13 Timer</b>	<b>35</b>
<b>14 Tips</b>	<b>36</b>
<b>15 Connections and interfaces</b>	<b>37</b>
15.1 RS 232 C serial interface	37
15.1.1 Configuration	38
15.1.2 Instruction syntax	38
15.1.3 Overview of the NAMUR-instructions	38
15.1.4 “Watchdog” -function, monitoring the serial data flow	41
15.2 Connections between RET control/t and external devices	41
15.2.1 PC 5.1 Adapter IKA - Control	41
15.2.2 PC 2.1 Cable	41
15.2.3 PC 1.5 Cable	41
15.2.4 PC 1.2 Adapter	42
15.2.5 AK 2.1 Cable	42
15.2.6 AK 2.2 Cable	42
15.3 Tips for using labworldsoft	42
<b>16 Faults</b>	<b>42</b>
16.1 Table error codes, troubleshooting, causes and solutions	43
<b>17 Maintenance and cleaning</b>	<b>46</b>
<b>18 Accessories</b>	<b>46</b>
<b>19 Technical data</b>	<b>47</b>

## 2 Safety instructions

EN

### *For your protection*

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Socket must be earthed (protective ground contact).
- **Caution - Magnetism!** Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).
- **Risk of burns!** Exercise caution when touching the housing parts and the heating plate. The heating plate can reach temperatures in excess of 340 °C. Pay attention to the residual heat after switching off.
- Ensure that the mains power supply cable does not touch the heating base plate.
- Wear your personal protective equipment in accordance with the hazard category of the media to be processed. Otherwise there is a risk from:
  - splashing and evaporation of liquids
  - ejection of parts
  - release of toxic or combustible gases.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Gradually increase the speed.
- Reduce the speed if
  - the medium splashes out of the vessel because the speed is too high
  - the appliance is not running smoothly
  - the container moves on the base plate.
- **Caution!** Never process and heat up any media that has a flash point higher than the adjusted safe temperature limit that has been set (50 to 350 °C).

The safe temperature limit must always be set to at least 25 °C lower than the fire point of the media used.

- Beware of hazards due to:
  - flammable materials
  - combustible media with a low boiling temperature
  - glass breakage
  - incorrect container size
  - overfilling of media
  - unsafe condition of container.
- The appliance may heat up when in use.
- The base plate can heat up due to the action of the drive magnets at high motor speeds, even if the heater is not operational.
- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA** if you have any questions.
- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Please observe the operating instructions for any accessories used.
- Ensure that the external temperature sensor (PT 1000, PT 100, ETS-D ...) is inserted in the media to a depth of at least 20 mm.
- The PT 1000 external temperature sensor must always be inserted in the media when connected.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- Always disconnect the plug before fitting accessories.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- The device will automatically restart in mode B following any interruption to the power supply.
- It may be possible for wear debris from rotating accessory parts


to reach the material being processed.

- When using PTFE-coated magnetic bars, the following has to be noted: *Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.*

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

#### For protection of the equipment

- The appliance may only be opened by experts.
- The voltage stated on the type plate must correspond to the mains voltage.
- Do not cover the device, even partially e.g. with metallic plates or film. This results in overheating.
- Ensure that the base plate is kept clean.
- Protect the appliance and accessories from bumps and impacts.
- Only use glass encapsulated temperature sensors for electrolysis procedures.
- Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 800 mm).

### 3 Correct use

#### • Use

- For mixing and/or heating liquids

#### • Range of use

- Laboratories
- Schools
- Pharmacies
- Universities

This device is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the correct use and the manufacturer's specifications or if the appliance is modified by third parties.

### 4 Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers:

a RET control/t, an operating instructions, a main cable and a short instruction.

## 5 Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). Insert the temperature measuring sensor (not included in the supply schedule) into the 5-polar DIN-bushing at the rear of the instrument. Put a magnetic bar into a vessel that is permeable for magnetic field lines, and fill the vessel with the medium to be processed. Put the vessel on the heating plate and immerse the temperature measuring sensor into the medium (min. 20 mm). See chapter 11.3 for the operation without external temperature measuring sensor. If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

## 6 Control elements




Fig. 1

## 7 Display

When the functions are switched off, the digital display always shows the set speed value.

In all three modes of operation only function "Heating" has been started, the display shifts to rated and real temperature.


In all three modes of operation only function "Stirring" has been started, the display shifts to rated and real speed.

If both functions have been started, the display of function „Heating“ always has preference. If in this case the speed is changed via rotary knob „Mot“ (not possible in **C**-mode of operation), the display reverses to speed display for a duration of 5 seconds.

The mode of the display can also be changed by pressing foil key „Mode“ once. After 5 seconds the display shifts back to temperature display.

If function „Heating“ is switched off, and if the heating plate temperature is higher than 50 °C, LED „Temp“ will be blinking green (20% ON ... 80% OFF).

The digital display changes constantly from rated temperature [°C]

to rated temperature

**HOT**

## 8 Switching on

The instrument is switched on by the mains switch on the left side of the instrument (see fig. 1).

A self test cycle for all the safety functions is run and a recognition cycle of the temperature sensors connected up is also run.

- The “Temp” and “Mot” LEDs and the LCD display segments flash for approximately 2 seconds.
- The software version number (e.g. “V000”) is displayed.
- The status of the intermittent mode (**On/OFF Intermit**) and the stirring bar breakaway monitoring (**On/OFF SAFE rpm**) functions are displayed.
- The message **TEST Probe** is displayed.
- The operating mode (**A, B** or **C**) is displayed.

- The temperature sensors connected up and identified are displayed:

**PT OFF** No external temperature sensor connected up.  
**PT2 OFF** PT100 (single)  
**PT2 On** 2x PT1000

- The safe temperature limit of the heating plate is displayed as **SAFE H xxx °C** (SV flashes).
- **SAFE H OK** is displayed.
- If external temperature sensors (PT 100 or PT 1000) are connected up, their safe temperature limit is displayed as **SAFE P xxx °C** (SV flashes).
- The message **SAFE P OK** is displayed.

## 9 Factory settings

The RET control/t is supplied with the following factory settings:

Operating mode	A	
Setpoint temperature	0	°C
SAFE H	350	°C
SAFE P	350	°C
Setpoint speed	50	rpm
Intermittent stirring mode	OFF	
Cycle time (at setpoint speed)	60	sec. (*)
Interruption (at 50 rpm)	10	sec. (*)
Stirring bar breakaway monitoring	OFF	
Error 5 response time	180	sec. (*)
PT 100 temperature sensor	0	K (*)
offset		
PT 1000 1 temperature sensor	0	K (*)
offset		
PT 1000 2 temperature sensor	0	K (*)
offset		
Timer		OFF

(\* = can only be modified via serial interface)

## 9.1 Resetting parameters to factory settings (RESET)

To reset the parameters to the factory settings, proceed as follows:

- Switch off the unit at the main ON/OFF switch.
- Hold the membrane buttons "Timer" and "Temp" pressed down and switch on the main ON/OFF switch. After approx. 10 seconds the display will begin to flash. Release the membrane buttons.
- The unit is now set on the factory settings.

## 10 Operation modes

The RET control/t can be operated in the following three modes of operation:

### 10.1 Operation mode A

The RET control/t is switched ex works to the operation mode **A**. In this mode, an operation is only possible when the **temperature measuring sensor** (PT 100 or PT 1000) is inserted. When the instrument is switched on („Power“), the functions "Heating" and "Stirring" are switched off.

The set rated values are stored and digitally displayed. They are taken over when the respective function is switched on. The set rated values can be changed. If there has been an interruption of the mains, **both functions have to be restarted**.

The safety circuit can be set or modified. When the power supply switch is switched on, **A** is displayed.

### 10.2 Operation mode B

In this mode of operation can be worked with or without **temperature measuring sensor** (PT 100 or PT 1000). When the device is switched on, the state of the functions "Heating" and "Stirring" before the last switching off (mains OFF) as well as the data set

last are taken over. The set rated values can be changed. The safety circuit can be set or modified.

The mode **B** of operation is indicated by symbol „**O**“.

When the power supply switch is switched on, **B** is displayed.

### 10.3 Operation mode C

In this mode of operation can be worked with or without **temperature measuring sensor** (PT 100 or PT 1000). When the instrument is switched on, the state of the functions "Heating" and "Stirring" before the last switching off (mains OFF), as well as the data set last are taken over. The set rated values from the operation modes **A** and **B** cannot be changed.

The safety circuit **cannot** be set or modified. When the power supply switch is switched on, **C** is displayed.

### 10.4 Operation mode X1

In the operation mode X1 the equipment behaves as in the operation mode **A** with the following exceptions:

- The allocation of the two temperature measuring sensors PT 1000 is firm.
- No Error 5, if the temperature measuring sensor PT 1000 is taken out of the reaction vessel.

### 10.5 Changing the operation modes

Changing the operation modes requires the following steps:

- Switch off the instrument by the mains switch.
- Keep foil key „Mot" and „Temp" pressed and switch on the device by mains switch, release the foil keys after approx. 1 second.

## 11 Function heating

The unit has a built-in 600 W heating plate with a stainless steel surface. Heating is controlled by a wear-proof, electronic TRIAC element.

The heating plate is kept at a constant temperature by a control circuit. The heating plate temperature can also be monitored from a separate, adjustable safety circuit. The two temperature sensors (PT 1000) for this purpose are built into the heating plate.

The heating plate safety circuit (SAFE H (hotplate)) can be set in the temperature range from 50° to 350 °C.

The unit can also be controlled from a single, external temperature sensor PT 100 or from twin temperature sensors PT 1000.

Any external temperature sensors must be plugged in before the unit is switched on so that the automatic self-test/sensor recognition cycle can take place.

The external temperature sensors (medium) are monitored by a separate, adjustable safety circuit. The medium safety circuit (SAFE P (probe)) can be set from 50 to 350 °C or to the maximum setting of the heating plate safety circuit.

Use the "Temp" knob to set the setpoint temperature required in the range from 0° to 340 °C (max.) The temperature set will be displayed on the upper right-hand side of the digital display. "SV" and "PROBE" will flash.

To switch on and off the function "Heating", use the foil key "Temp".

When the function "Heating" is switched on, the green LED "Temp" will light up and the digital display will switch to the next display page. The actual temperature is displayed.

While the plate is heating up, the green LED "Temp" will change colour from green to orange (orange = power supply to heating plate).


### 11.1 Medium temperature control by temperature sensor PT 100

When the unit is switched on, the message "PT2 OFF" is displayed.

The "PROBE" symbol is shown on the digital display to indicate that the sensor PT 100 is operating. The setpoint and actual temperatures displayed are for the medium.

The heating plate temperature is limited by the setting of the safety temperature circuit (SAFE H).

The medium temperature is controlled and regulated by a fuzzy logic controller. This ensures that the medium, sensed by the sensor PT 100, is heated as quickly as possible without exceeding to the set temperature. The fuzzy logic controller automatically adjusts to the various heating capacities of different vessels and their content. This system ensures optimum temperature control without temperature drift or fluctuations.

### 11.2 Medium temperature control by twin temperature sensors PT 1000

When the unit is switched on, the message "PT2 OFF" is displayed.

The "PROBE" symbol is shown on the digital display to indicate that two sensors PT 1000 are operating.

Heat transfer medium and reaction medium (two vessels)

There is a temperature sensor in the heat transfer medium and a temperature sensor in the reaction medium.

Automatic recognition of the heat transfer and reaction medium sensors takes place.


The setpoint temperature displayed is for the reaction medium temperature. The actual temperature of the reaction medium is shown in large figures to the left of the setpoint temperature.

The heating plate temperature and the heat transfer medium temperature are limited by the heating plate safety temperature (SAFE H) setting. The reaction medium temperature is limited by the medium safety temperature (SAFE P) setting.

The reaction medium temperature is controlled and regulated by a fuzzy logic controller. The controller is optimised by the heat transfer medium temperature. This ensures that the reaction medium is heated as quickly as possible without exceeding the set temperature. The fuzzy logic controller automatically adjusts to various heating capacities of different vessels and their content. This system ensures optimum temperature control without temperature drift or fluctuations.

### 11.3 Medium temperature control without external temperature sensors

This is not possible in operation mode **A**. Medium temperature control using external temperature sensors is more precise than standard heating plate temperature control. The operation without PT 100 is digitally shown by the symbol "PROBE" blinking. The display rated and real temperature refer to the temperature of the heating plate.

### 11.4 Residual heat warning (HOT)

In order to prevent the risk of burns from a hot heating plate, the RET control/t has a residual heat (HOT) warning function. When the heating function is switched off and the heating plate temperature is still above 50 °C, the green LED "Temp" will flash. The digital display will switch backwards and forwards between "Setpoint temperature+ °C" and "Setpoint temperature + HOT" as a warning that the plate is still hot.

If the main ON/OFF switch on the unit is switched off, the display will remain active and the message "—" will flash in the area of the actual value display area. In this state the stirring function can be switched on and the setpoint speed can also be modified.

If the heating plate temperature drops to below 50 °C, the unit will automatically switch off. The display will be switched off and the unit electronics will be powered down by a relay. If you have to switch the heating plate off suddenly when the heating plate is still hot, you must switch the main ON/OFF switch on and then off again rapidly.

The residual heat warning function will not operate in the event of a mains power failure or if the power plug is pulled out of the socket.

## 11.5 Changing the temperature unit of measurement from °C to °F

To toggle the unit of measurement for temperature readings from °C to °F and back again, proceed as follows:

- Switch off the unit at the main switch ON/OFF. (The residual heat display should no longer be enabled, see residual heat Warning (HOT)).
- Hold the foil keys "Mode" and "Temp" pressed down and switch on the main switch ON/OFF. After approx. 5 seconds the display will begin to flash. Release the foil keys.
- Temperatures will now be displayed in the unit of measurement selected. Temperatures in the memory will be converted accordingly.

## 11.6 Setting the safety temperature

Setting and/or changing the safety temperature is only possible in the operation mode **A** or **B**.

The safety circuit is set as follows:

Switch on the device by mains switch on the left side of the instrument.

As soon as the display shifts to SAFE H (safety temperature limiter of the heating plate) (SV is blinking), press the foil key "Mode" (SV is blinking no longer), keep this foil key pressed and set the requested safety temperature with the rotary knob TEMP. The set safety temperature (e.g. 100 °C) is displayed. After releasing the foil key "Mode", SAFE H OK is displayed for about 2 seconds, and the safety temperature is stored this way.

When the display changes to show SAFE P (safety temperature limit for external sensors) (SV flashes), press the foil key "Mode". The SV will stop flashing. Hold down the foil key "Mode" and turn the knob TEMP to set the safety limit temperature required. The safety temperature (e.g. 60 °C) will be displayed. Release the foil key "Mode". After approx. 2 seconds, the message SAFE OK will be displayed to indicate that the safety temperature has now been stored in the memory.

## 12 Function stirring

The motor speed is set on the rotary knob "Mot", (0 and/or 50 to 1200 rpm in steps of 10 rpm each). "SV" is blinking.

The function is started and stopped with foil key "Mot".

If the function "Stirring" is started, the motor runs smoothly to the set speed, and the LED "Mot" beside the imprint "Mot rpm" shines green.

The display changes to real speed display.

### 12.1 Intermittent mode

If Intermittent mode is active, "On IntErmit" will be displayed when the device is switched on.

The intermittent mode interrupts the stirring cycle at regular intervals. For example, if the function "Stirring" is switched on with the setpoint speed set at 1 rpm the stirring cycle will be interrupted for approx. 10 seconds after a cycle time of 60 seconds. The cycle time and the length of the interruption can be modified via the serial interface (see section 15.1).

**The stirring bar breakaway monitoring function cannot be enabled when the intermittent mode is enabled.**

The intermittent mode can also be used with the special small stirring bars. These can be used when the stirring bar breakaway monitoring function is not enabled and will continue to stir the medium even when the stirring bar breaks away. These stirring bars are picked up once again when the speed drops cyclically.

## 12.2 Stirring bar breakaway monitoring

When the unit is switched on, the stirring bar breakaway monitoring function is displayed (On/OFF SAFE rpm).

Stirring bar breakaway can occur with any magnetic stirrer. Breakaway can occur when stirring at high speeds or when large stirring bars are being driven at speeds close to the magnetically transmissible maximum torque and there are variations in medium viscosity or controller fluctuations.

On the RET control/t you can select the stirring bar breakaway monitoring function which operates in the setpoint speed range between 300 and 1200 rpm.

When stirring bar breakaway is detected, stirring is stopped momentarily to enable the stirring bar to be picked up. After this, stirring starts again and ramps up to the setpoint speed (< flashes). If a further stirring bar breakaway occurs within the next 3 minutes of stirring at a constant setpoint speed, the stirring bar will be picked up again and stirring will restart but this time at the setpoint speed reduced by 100 rpm. The red LED "Temp" flashes and the buzzer will sound for several seconds. After the fourth stirring bar breakaway – with less than 3 minutes normal stirring – the error message ER 41 will be displayed (the motor will continue to run and the heating is switched off) the red LED "Temp" will flash and the buzzer will sound. When further stirring bar breakaways occur, the speed will be reduced each time in 100 rpm steps until the minimum of 200 rpm is reached. The speed is no longer displayed. In this case you should select new test conditions.

**IMPORTANT: You must calibrate and test the stirring bar breakaway monitoring function for your particular configuration taking into account the stirring bars, speed, vessel and medium you are using.**

The reliability of stirring bar breakaway monitoring increases as the transmitted torque increases.

In a glass beaker filled with water, the function will detect the breakaway of stirring bars longer than 50 mm (d=8 mm) in the speed range from 300 to 1200 rpm. Increasing the stirring bar length, the stirring bar surface area (e.g. TRIKA coated bars) or increasing the speed will increase the torque transmitted.

**The stirring bar breakaway monitoring function will not operate if the stirring bar breaks away while the stirring is still ramping up to the setpoint speed.**

## 12.3 Setting the intermittent mode and the stirring bar breakaway monitoring function

- Switch off the unit at the main switch ON/OFF.
- Hold the foil keys "Mot" and "Mode" pressed down and switch on the main switch ON/OFF. After approx. 5 seconds release the foil keys.
- After the display self-test (when all the LCD segments flash), the message On or OFF IntErmitt will be displayed.
- Press the foil key "Mot" within the next 5 seconds to switch the intermittent mode ON or OFF as required.
- If you disable the intermittent mode by selecting OFF, the message On or OFF Safe rpm will be displayed.
- To switch the stirring bar monitoring function ON or OFF, press the foil key "Mot" within the next five seconds.

**IMPORTANT: These changes will be stored to memory when the self-test (SAFE H or SAFE P OK is displayed) has been completed.**

## 13 Timer

Set the function "Heating" by the timer.

- Switching on the timer  
Activate the timer by pressing the "Timer" membrane key; the "SV" symbol will flash.
- Setting the timer  
Set the desired time (0:00 - 9:59 h) by turning the left knob minute by minute.  
Caution! When the time is set to 0:00 the timer function is switched off!
- Starting the timer  
Start the timer by switching on the "Heating" function; the "⌘" symbol will flash.

Caution! Always make sure the "Heating" function is switched off before starting the timer. Then switch on the timer followed by the "Heating" function.

- Displaying the timer

Press the "Mode" membrane key to display the remaining time. The "Temperature", "Speed" and "Timer" values will be displayed. If you switch on the "Stirring" function when the timer is activated, the timer will not restart.

To switch on the "Stirring" function, press the "Mot" membrane key.

- Changing the timer

Activate the timer by pressing the "Timer" membrane key.

Set the desired time (0:00 - 9:59 h) by turning the left knob minute by minute.

Confirm by switching the "Heating" function off and back on again.

- Stopping the timer

To stop the timer during countdown, switch the "Heating" functions off. You can then switch the heating and/or stirring function(s) back on and the DIGITAL-T will function in continuous mode.

Once the countdown has finished, an audible signal will be given, the "Heating" function will be switched off, the "Stirring" function is still active. The display will remain lit up for as long as the power supply is on.

- Timer in modes B and C

If the appliance was switched off before the countdown had finished, the timer will run for the remaining time the next time the appliance is switched on.

If the appliance was switched off after the timer had finished, the set values will be retained, however the "Heating" and "Stirring" functions will be switched off.

Caution! To switch off the timer in C mode, switch off the "Heating" function. Switch the "Heating" function back on without activating the timer function or set the timer to 0:00.

- Timer in mode X1

In mode X1 the DIGITAL-T functions the same way as in mode A with the following differences:

The allocation of the two PT 1000 temperature sensors (heat carrier and reaction medium) is fixed.


No error 5 when the PT 1000 temperature sensor is removed from the reaction vessel (the temperature of the sensor continues to be displayed). The heat carrier sensor continues to be monitored.

## 14 Tips

- Ensure that the medium is well mixed (high speeds and large stirring bars ensure good mixing). A medium which is well mixed will ensure uniform temperature distribution and better temperature control.
- The temperature sensors (PT 100, PT 1000) should be immersed to the maximum possible depth. The minimum immersion depth is 20 mm; do not go below this immersion depth.
- When using metal vessels, do not place the temperature sensors on the bottom of the vessel. Placing sensors on the vessel bottom can cause excessively high temperatures to be measured especially in media which have poor conductivity. The tip of the measuring sensor must be at least 5 mm from the vessel bottom; a distance of 10 mm is ideal.
- With highly viscous reaction media, the distance between the heating plate and the stirring bars must be kept as small as possible (to reduce the risk of stirring bar breakaway).

## 15 Connections and interfaces

### 15.1 Serial RS 232 C interface

#### 15.1.1 Configuration

The unit can be controlled from an external PC (using the dedicated software labworldsoft) via the RS 232 C serial interface fitted to the unit.

**To increase safety when controlling the RET control/t from a PC, enable the watchdog function for monitoring the continuous serial data flow (see section 15.1.4).**

- The functions of the interface lines between laboratory instrument and automation system are selected from the specified signals of the EIA-standard RS 232 C, corresponding with DIN 66 020 Part 1. The allotment of the bushing can be taken from illustration.
- For the electrical properties of the interface lines and for the allotment of the signal status, standard RS 232 C, corresponding with DIN 66 259 Part 1 applies.
- Transmission method: Asynchronous signal transmission in start-stop-operation.
- Mode of transmission: Fully Duplex.
- Character format: Character imaging acc. to data format DIN 66 022 for start-stop-operation. 1 start bit; 7 character bits; 1 parity bit (straight (even)); 1 stop bit.
- Transmission speed: 9600 bit/s
- Data flow control:
  - Hardware handshake RTS/CTS
  - RTS: (pin 7) LOW (positive tension) / computer may send
  - RTS: (pin 7) HIGH (negative tension) / computer may not send
  - CTS: (pin 8) LOW (positive tension) / computer recipient
  - CTS: (pin 8) HIGH (negative tension) / computer not recipient
- Access method: Data communication from laboratory instrument to computer is only possible on demand of the computer.

### 15.1.2 Instruction syntax

Here applies the following:

- The instructions are generally sent from the processor (master) to the laboratory instrument (slave).
- The laboratory instrument exclusively sends on demand of the processor. Even error codes cannot be spontaneously communicated from the laboratory instrument to the processor (automatic system).
- The instructions are transmitted in capital letters.
- Instructions and parameters as well as subsequent parameters are separated by at least one blank (Code: hex 0x20).
- Each individual instruction including parameters and data as well as each reply are terminated with CR LF (Code: hex 0x0D and 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a floating point number is the point (Code: hex 0x2E).

The above statements largely correspond with the recommendations of the NAMUR-Association. (NAMUR-recommendations for the design of electric plug connections for the analog and digital signal transmission to laboratory - MSR individual units. Rev. 1.1).

### 15.1.3 Overview of the NAMUR-instructions

Abbreviations:

X,y	=	Numbering parameter (integer number)
m	=	Variable value, integer
n	=	Value of variable, floating point number
X = 1	=	PT 100 or PT 1000 medium temperature
X = 2	=	Heating plate temperature
X = 3	=	Heating plate - safety temperature
X = 4	=	Speed
X = 7	=	PT 1000 heat transfer medium temperature
X = 8	=	PT temperature sensor – safety temperature
X = 50	=	PT 100 temperature sensor offset in K ( -3.0 <= n <=+3.0 )
X = 52	=	PT 1000 1st temperature sensor offset in K ( -3.0 <= n <=+3.0 )
X = 53	=	PT 1000.2nd temperature sensor offset in K ( -3.0 <= n <=+3.0 )
X = 54	=	Error 5 response time in seconds ( 180 <= n <= 1200 )
X = 55	=	Intermittent mode cycle time in seconds ( 10 <= n <= 600 )
X = 56	=	Intermittent mode, duration of interruption in seconds ( 5 <= n <= 60 )

NAMUR instructions	Function	Display (additional)
IN_NAME	Input description name	
IN_PV_X      X=1;2;3;4;5; 7; 8;	Reading the real value	
IN_SOFTWARE	Input software ID number, date and version	
IN_SP_X      X=1;2;3;4;7;8; 12;42;50;52;53; 54;55;56;	Reading the set rated value	
IN_TYPE	Input laboratory unit ID	
OUT_NAME name	Output description name. (Max. 6 characters, default: IKAMAG)	
OUT_SP_12@n	Setting the WD safety temperature with the echo of the set value	
OUT_SP_42@n	Setting the WD safety speed with the echo of the set value	
OUT_SP_X n    X=1;2;4;7;50; 52;53;54;55;56	Setting the rated value to n	
OUT_WD1@m	Watchdog mode 1: When a WD1 event occurs, the heating and stirring functions are shutdown and message ER 2 is displayed. Set the watchdog time to m (20...1200) seconds, with echo of the watchdog time. This instruction starts the watchdog function and must be sent within the set watchdog time.	
OUT_WD2@m	Watchdog mode 2: When a WD2 event occurs, the speed setpoint will be set to the WD safety setpoint speed and the temperature setpoint will be set to the WD safety setpoint temperature. The WD warning is displayed. The WD2 event can be reset with OUT_WD2@0 - resetting also blocks the watchdog function. Set the watchdog time to m (20...1200) seconds, with echo of the watchdog time. This command starts the watchdog function and must be sent within the set watchdog time.	
RESET	Switching off the instrument function.	

RMP_CONT_X X=1;4	Continuation of ramp function (After prior RMP_PAUSE_X).	
RMP_IN_X X=1;4	Reading the real segment number of ramp. With ramp not started: 0	
RMP_INX_y X=1;4	Reading the accumulated value and the ramp segment duration of ramp segment y.	
RMP_LOOP_RESET_X X=1;4	Ending of ramp loop.	
RMP_LOOP_SET_X X=1;4	To work off the ramps in one loop.	
RMP_OUT_X_yn X=1;4 hh:mm:ss	Setting the accumulated value (n) and the ramp segment duration (hh:mm:ss) for ramp segment y.	
RMP_PAUSE_X X=1;4	Stopping the ramp function. Freezing of real rated value and real ramp segment time.	
RMP_RESET_X	Switching off ramp functions and deleting of all set ramp segments	
RMP_START_X X=1;4	Starting the ramp function, beginning with ramp function No.1 (Only possible after prior START_X. After RMP_STOP_X START_X is not necessary.	
RMP_STOP_X X=1;4	Switching off ramp function. Rated value = 0. (Ramp is maintained, that means, ramp can be restarted with RMP_STARTX_X).	
START_X X=1;2;4;5;7	Starting the instrument's (remote) function.	REMOTE
START_X X=16 (old function)	Activation of the "watchdog" function. If no rated value is entered given in the next 60 seconds, the heating is switched off.	
STATUS_X	<p>Display of status</p> <p>1S: mode of operation A</p> <p>2S: mode of operation B</p> <p>3S: mode of operation C</p> <p>S0: manual operation without fault</p> <p>S1: Automatic operation Start (without fault)</p> <p>S2: Automatic operation Stop (without fault)</p> <p>&lt;0: error code: (-1)</p> <ul style="list-style-type: none"> <li>- 1: error 1</li> <li>- ... (see table)</li> <li>-31: error 31</li> <li>-83: wrong parity</li> <li>-84: unknown instruction</li> <li>-85: wrong instruction sequence</li> <li>-86: invalid rated value</li> <li>-87: not sufficient storage space</li> </ul>	
STOP_X X=1;2;4;5;7	Switching off the instrument - (remote) function. Variables set with OUT_SP_X are maintained. Contains the instruction RMP_STOP.	REMOTE

### 15.1.4 "Watchdog" function, monitoring the serial data flow

The following applies to situations where the watchdog function is enabled (see Namur instructions). If no new transmissions of these commands from the PC take place within the preset watchdog time, the heating and stirring functions will be shutdown according to the watchdog mode selected or will be controlled using the preset setpoints. An operating system crash, a PC power failure or a fault in the connecting cable to the RET control/t can cause an interruption in data transmission.

#### "Watchdog" Mode 1

If an interruption in data transmission occurs which is longer than the preset watchdog time, the heating and stirring functions will be shutdown and the error message ER 2 will be displayed.

#### "Watchdog" Mode 2


If an interruption in data transmission occurs which is longer than the preset watchdog time, the speed setpoint value will be set to the WD safety speed setpoint and the temperature setpoint will be set to the WD safety temperature setpoint. The WD warning message will be displayed.

## 15.2 Connections between DIGITAL-T and external devices

The RET control/t unit can be connected to other external devices using the following IKA adapters and adapter cables.


### 15.2.1 PC 5.1 Adapter IKA - Control

An adapter cable is available as an accessory part. It splits the analog and serial signals. The analog output signals are allotted to a 7-polar coupling according to the NAMUR-recommendation; the serial signals to a 9-polar Sub-D-bushing plug (RS 232 C).

### 15.2.2 PC 2.1 Cable

Cable PC 2.1 for connection to the 9-pin socket on the computer.

### 15.2.3 PC 1.5 Cable

Cable PC 1.5 for connection to the 15-pin socket on the computer (25-pin plug).

### 15.2.4 PC 1.2 Adapter

Adapter PC 1.2 for connecting the 9-pin socket with an 8-way serial interface (25-pin plug).

### 15.2.5 AK 2.1 Cabel

Cable AK 2.1 for connecting the 7-pin socket of the PC 5.1 adapter to a printer (4 mm banana plug).

### 15.2.6 AK 2.2 Cabel

Cable AK 2.2 for connecting the 15-pin socket to a printer (4 mm banana plug).

## 15.3 Tips for using labworldsoft

- The safety temperatures (SAFE H, SAFE P) must be set to match the test conditions set on the RET control/t also when you are using labworldsoft.
- If the setpoint settings in labworldsoft exceed the setpoint value range or the safety temperature set on the RET control/t, the labworldsoft settings will not be accepted. The speed setpoint must not exceed 1200 rpm.
- If the ramp function is used to increase gradually the medium setpoint temperature, the system variation and the control hysteresis can be greater than they would be if a setpoint temperature is entered directly.
- Heating plate temperature control in critical processes should only be used by experienced users. Critical processes in this case refers to the testing of large volumes under poor control conditions with a single medium temperature setpoint setting.
- When a watchdog mode is enabled, the watchdog time must be set at 20 times the sensing time. This is to prevent process delays, caused by parallel programs (e.g. Excel) or devices (e.g. NPort-Server), from triggering an Error 2 (watchdog time timeout).

## 16 Faults

A unit fault is indicated by the red flashing LED "Temp" (flashing in one second intervals) and the display of an error code number "Er n" on the display (see the error code table). The buzzer also sounds.

To switch off the warning buzzer, press the foil key "Mode".

The heating will be powered down (exception: WD – see table).

Following a fault in the heating circuit, the stirring function will continue to operate at the speed set before the fault took place (exception: WD - see table).

To reset a fault, switch off the unit at the main switch ON/OFF for a few seconds. The unit should remain off long enough for the temperature to settle (the residual heat display must be OFF, see the residual heat display).

When the unit is switched on again, the safety circuit test is run within 20 seconds. After this the other test/monitoring functions are enabled.

**If the fault will not clear, you must run a RESET cycle (see section 9.1 "Resetting parameters" to factory settings (RESET)).**

**If a RESET cycle does not clear the fault, take the device to your technical service center.**

## 16.1 Troubleshooting, Error messages, causes and solutions

Error code	Description	Cause	Effect	Solution
Er 1	In operating mode A; temperature sensor not plugged in	<ul style="list-style-type: none"> <li>When the unit was switched on in the mode A, the temperature sensor was not plugged in</li> </ul>	The heating remains switched off	<ul style="list-style-type: none"> <li>Switch off the unit, plug in the temperature sensor, switch on the unit</li> <li>Switch over the operating mode</li> </ul>
Er 2	In remote operation (PC), with watchdog funktion 1 enabled: No communication between PC and RET control/t	<ul style="list-style-type: none"> <li>PC does not send any data during the watchdog time.</li> <li>The connection/cable to the PC is broken.</li> </ul>	Heating OFF Motor OFF	<ul style="list-style-type: none"> <li>Change the watchdog time.</li> <li>Send data from the PC within the watchdog time (OUT_WDx@m).</li> <li>Check the connector cable and plug.</li> </ul>
Er 3	Temperature inside unit is higher than 80 °C	<ul style="list-style-type: none"> <li>Build up of heat between the heating plate and the casing.</li> <li>Permitted ambient temperature has been exceeded.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>Switch off the unit. Allow it to cool down and then switch on again.</li> <li>Change test process.</li> <li>Restore ambient temperature to within the permitted maximum.</li> </ul>
Er 4	Difference between setpoint and actual speeds is greater than 300 rpm.	<ul style="list-style-type: none"> <li>Motor seized or overloaded.</li> </ul>	Heating OFF Motor OFF	<ul style="list-style-type: none"> <li>Reduce the torque load, e.g. use smaller magnetic bars.</li> <li>Reduce the setpoint speed.</li> </ul>
Er 5	No temperature increase measured by temperature sensor where temperature difference remains the same (after approx. 3 minutes).	<ul style="list-style-type: none"> <li>Sensor not in medium.</li> <li>Volume of heat transfer medium too large.</li> <li>Medium heat transfer rating is too low.</li> <li>Heat temperature rating of the vessel is too low.</li> <li>With indirect heating the total heat transfer resistance is too high.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>Immerse the sensor in the medium.</li> <li>Reduce medium volume.</li> <li>Use a heat transfer oil with a better heat transfer rating.</li> <li>Use a metal vessel in place of a glass vessel.</li> </ul>
Er 8	The calibration value of the temperature sensor is outside the limit value.	<ul style="list-style-type: none"> <li>Fault in calibration procedure.</li> <li>Value was incorrectly stored to memory.</li> <li>EEPROM switching error.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>Repeat the calibration procedure.</li> </ul>
Er 9	Fault storing the safety circuit setpoint temperatures.	<ul style="list-style-type: none"> <li>EEPROM fault</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 11	Temperature sensor unplugged.	<ul style="list-style-type: none"> <li>Medium temperature sensor 1 or 2 (PT 100 or PT 1000) was unplugged during an operating cycle.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>Only unplug the temperature sensor when the unit is switched off.</li> </ul>

<b>Error code</b>	<b>Description</b>	<b>Cause</b>	<b>Effect</b>	<b>Solution</b>
Er 12	Temperature sensor plugged in.	<ul style="list-style-type: none"> <li>• Medium temperature sensor 1 or 2 (PT 100 or PT 1000) was unplugged during an operating cycle.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Only unplug the temperature sensor when the unit is switched off.</li> </ul>
Er 13	SAFE TEMP P (Probe) adjustable safety circuit fault.	<ul style="list-style-type: none"> <li>• Difference between the setpoint and actual values of the adjustable safety circuit for the minimum temperature monitoring.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Switch ON and set the SAFE TEMP P to another value. If this clears the fault when the unit is switched on again, you can return to using the previous value.</li> </ul>
Er 14	Temperature sensor 1, short-circuit	<ul style="list-style-type: none"> <li>• Short-circuit at temperature sensor plug.</li> <li>• Short-circuit in connecting cable or on temperature sensor.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Check the plug.</li> <li>• Substitute the temperature sensor.</li> </ul>
Er 15	Temperature sensor 2, short-circuit	<ul style="list-style-type: none"> <li>• Short-circuit at temperature sensor plug.</li> <li>• Short-circuit in connecting cable or on temperature sensor.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Check the plug.</li> <li>• Substitute the temperature sensor.</li> </ul>
Er 16	Temperature sensor 1 has exceeded the SAFE TEMP P (Probe)	<ul style="list-style-type: none"> <li>• SAFE TEMP P (Probe) is set at a temperature lower than the actual temperature on the temperature sensor 1.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Allow the medium to cool.</li> <li>• Set the SAFE TEMP P (Probe) at a higher temperature</li> </ul>
Er 17	Temperature sensor 2 has exceeded the SAFE TEMP P (Probe).	<ul style="list-style-type: none"> <li>• SAFE TEMP P (Probe) is set at a temperature lower than the actual temperature on the temperature sensor 2.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Allow the medium to cool.</li> <li>• Set the SAFE TEMP P (Probe) to a higher temperature</li> </ul>
Er 21	Fault during heating plate safety test.	<ul style="list-style-type: none"> <li>• Safety relay does not open.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 22	Fault during heating plate safety test.	<ul style="list-style-type: none"> <li>• S_CHECK cannot generate a higher H_S_TEMP.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 23	Fault in SAFE TEMP H (Hotplate) adjustable safety circuit	<ul style="list-style-type: none"> <li>• Difference between the setpoint and actual values of the adjustable safety circuit for the heating plate temperature monitoring</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Switch ON in operating mode A and set the SAFE TEMP P to another value. If this clears the fault when the unit is switched on again, you can return to using the previous value.</li> </ul>

Error code	Description	Cause	Effect	Solution
Er 24	Heating plate temperature has exceeded the SAFE TEMP H (Hotplate)-setting	<ul style="list-style-type: none"> <li>• SAFE TEMP H (Probe) is set a temperature lower than the heating plate actual temperature.</li> <li>• Breakage of the heating plate temperature control sensor</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Allow the heating plate to cool.</li> <li>• Set SAFE TEMP H (Hotplate) to a higher temperature.</li> </ul>
Er 25	Heating switching element monitoring	<ul style="list-style-type: none"> <li>• There is a short-circuit in a switching element (Triac) of the heating control circuit</li> <li>• Safety relay has shutdown the heating circuit.</li> <li>• The heating or the connection is broken</li> <li>• The heating plate safety temperature sensor is broken.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Switch off the unit. Allow it to cool down and then switch on again.</li> </ul>
Er 26	Difference between the internal heating plate control and safety temperature sensors is too large	<ul style="list-style-type: none"> <li>• Irregular temperature distribution on the heating plate due to point heat dissipation</li> <li>• Fault in the control or safety temperature sensors.</li> </ul>	Heating OFF	<ul style="list-style-type: none"> <li>• Switching off the unit. Allow in to cool down and then switch on again.</li> <li>• When using metal blocks on the heating plate surface, ensure regular heat dissipation</li> </ul>
Er 31	Fault in the heating switch element (Triac)		Heating OFF	
Er 41	Stirring bar breakaway	<ul style="list-style-type: none"> <li>• Stirring bar has broken away on four consecutive occasions without any changes to the setpoint value.</li> </ul>	The setpoint speed is reduced to 300 rpm (min.) Heating OFF	<ul style="list-style-type: none"> <li>• Select a lower speed setpoint.</li> <li>• Kleineres Rührstäbchen wählen</li> </ul>
WD (Watch Dog)	In remote operation (PC) with watchdog function 2 enabled: No communication between PC and RET control/t	<ul style="list-style-type: none"> <li>• PC does not send OUT_WDx@m data during the watchdog time</li> <li>• The connection/cable to the PC is broken</li> </ul>	The setpoint temperature is set to the WD safety temperature The setpoint speed is set to the WD safety speed	<ul style="list-style-type: none"> <li>• Change the watchdog time.</li> <li>• Send data from the PC within the watchdog time (OUT_WD@m).</li> <li>• Check the connector cable and plug.</li> </ul>

## 17 Maintenance and cleaning

The RET control/t is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

Use only cleansing agents which have been approved by IKA to clean IKA devices.

To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA.

Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

### Spare parts order

When ordering spare parts, please give:

- machine type,
- manufacturing number, see type plate,
- item and designation of the spare part, see [www.ika.net](http://www.ika.net), spare parts diagram and spare parts list.

**Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.**

Send in the case of service the equipment back in the original packaging zurück. Storage packing is not sufficient for the back dispatch. Use additionally a suitable transportation packing.

## 18 Accessories

- PT 100.50 PT 100 Temperature measuring sensor, stainless steel, DIN EN 60751 Cl. A
- PT 100.51 PT 100 Temperature measuring sensor, glass-coated, DIN EN 60751 Cl. A
- PT 1000.50 Two-way PT 1000 Temperature measuring sensor, stainless steel, DIN EN 60751 Cl. A
- PT 1000.51 Two-way PT 1000 Temperature measuring sensor, class-coated, DIN EN 60751 Cl. A
- H 15 Bath top stainless steel 1,5 l
- H 16V Support rod
- H 36 Holding rod
- H 44 Cross sleeve
- H 28 Bath top stainless steel 1 l
- H 29 Oil bath top aluminium 1 l
- H 30 Oil bath top aluminium 1,5 l
- RS 1 Set of stirring bars
- RSE PTFE-stirring bar remover
- PC 5.1 Adapter
- PC 1.2 Adapter
- PC 1.5 Cabel
- PC 2.1 Cabel
- AK 2.1 Cabel
- AK 2.2 Cabel

- labworldsoft

### Recommended IKA-stirring bars

- TRIKA-stirring bar PTFE, triangular, lengths 25, 40 and 55 mm
- IKAFLON-stirring bar PTFE
  - ø 6mm lengths 10 and 15 mm
  - ø 8mm lengths 20, 25, 30, 40, 50 and 60 mm

## 19 Technical data

Design voltage	VAC	220...230 ±10%
or	VAC	115±10%
or	VAC	100±10%
Design frequency	Hz	50/60
Display		digital
Display dimensions	mm	47 x 21
Timer	h	0:00 - 9:59
<u>Stirring</u>		
Motor		EC-motor with ball bearings, without brushes
Max. power input of motor	W	12
Max. power output of motor	W	5
Speed range	rpm	0; 50 at 1200
Setpoint speed setting resolution	rpm	10
Actual speed display resolution	rpm	1
Speed variation	rpm	±30
(no load, 115/230 V, at 1200 rpm, RT 25 °C)		
Stirred quantity (H <sub>2</sub> O)	l	20
<u>Heating</u>		
Heating plate material		Stainless steel
Ø of the heating plate	mm	135
Heating power (115/230 V)	W	600 ±10%
(100 V)	W	455 ±10%
Temperature range	°C	RT ...340
Setpoint temperature setting resolution	K	0,5 <100 °C
	K	1 >=100 °C
Actual temperature solution		
- medium temperature	K	0,1 <100 °C
	K	1 >=100 °C
- heat transfer temperature	K	1
Heating plate adjustable safety circuit	°C	50 - 350
Medium plate adjustable safety circuit	°C	50 - 350
Safety circuit setting resolution	K	1

Temperature sensor in medium		PT 100 or 2xPT 1000, DIN EN 60751 Kl. A
Temperature sensor - variation (PT 100 oder 2xPT 1000, DIN EN 60751 Kl. A)	K	$<= \pm(0,15 + 0,002^{\circ} T )$
Temperature measuring range variation	K	$<= \pm 0,3 < 100 \text{ }^{\circ}\text{C}$
	K	$<= \pm 0,5 > 100 \text{ }^{\circ}\text{C}$
Temperature measuring range + sensor variation	K	$<= \pm 0,45 \text{ at } 0 \text{ }^{\circ}\text{C}$
	K	$<= \pm 0,85 \text{ at } 100 \text{ }^{\circ}\text{C}$
Heating plate temperature variation	K	$\pm 3$ without vessel, heating plate centre at $100 \text{ }^{\circ}\text{C}$
Control hysteresis with temperature sensor (500 ml water, 600 ml glass beaker, stirring bar 30, 800 rpm, PT100.50, $50 \text{ }^{\circ}\text{C}$ )	K	$\pm 0,2$ (FUZZY-optimized)
Heating plate control hysteresis	K	$\pm 3$ without vessel, heating plate centre at $100 \text{ }^{\circ}\text{C}$
Heating speed for 1 litre water in H15 ( $25 \text{ }^{\circ}\text{C}$ at $99 \text{ }^{\circ}\text{C}$ , 600 W)	min	ca. 16
Analog output: 0...10 V		
Temperature		1V/100 K
Speed		1V/1000 rpm
Analog output voltage variation	V	$\pm 0,05$
Permitted ambient temperature	$^{\circ}\text{C}$	+5 at +40
Permitted relative humidity	%	80
Permitted duration of operation	%	100
Protection class acc. DIN EN 60529		IP 42
Fuses		F1/F2 6,3AT/250V Id-Nr 1408000 F3 1 A (Trafo secondary 20 Vac) F4 0,75 A (Trafo secondary 10 Vac)
Contamination level	2	
Overvoltage category		II
Protection class		1 (protective earth)
Operation at a terrestrial altitude	m	max. 2000 above sea level
Dimensions	mm	160 x 280 x 100
Weight	kg	2,8

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GmbH & Co.KG**
**IKA® Works Inc.**