IKAMAG® RO5power / RO10power / RO15power

IKAMAG® RT5power / RT10power / RT15power

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OPERATING INSTRUCTIONS EN 7

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DE

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CE-DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EEC; 92/31/EEC und 73/23/EEC and conforms with the standards or standardized documents EN 61010-1 (IEC/CEI 61010-1) and EN 61326-1 (IEC 61326-1).

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Warranty

You have purchased an laboratory machine which meets the highest engineering and quality standards.

In accordance with guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

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Contents

	Page
Warranty	2
Safety instructions	7
Correct use	7
Unpacking	8
Commissioning	8
Function Stirring	8
Function Heating	8
Maintenance and cleaning	9
Associated standards and regulations	9
Accessories	9
Technical data	9
List of spare parts RO / RT power	10
Spare parts diagram RO / RT power	15

Safety instructions

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening.

When using PTFE-coated magnetic bars, please pay attention to the following: *Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodic system at temperatures above 300 - 400°C. Only elementary fluorine, chlorine trifluoride and alkaline metals do attack PTFE,*

halogen hydrocarbons have a reversibly swelling effect.

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Attention: Magnetism! Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac, cardiac pacemakers...).

If the power cable is replaced, it must be by one of the same quality. Please put the supplied transparent foil on the supporting surface. The foil prevents that the glass vessels don't slip to the side.

RT: Caution: Never heat up any medium with ignition point higher than 180°C.

The instrument may not be operated in rooms with explosion hazard.

Furthermore, it has to be placed upon a fire-proof and/or not burnable surface.

RT: Please make sure that the mains cable does not contact the heating plate!

RT: Exercise caution when touching the housing parts and the heating plate!

Risk of burns!

The heating plate can reach temperatures in excess of 125°C.

RT: Please consider possible endangering due to release of toxic or flammable gas caused by heated medium.

Correct use

IKAMAG RO 5 power, IKAMAG RO 10 power and IKAMAG RO 15 power are laboratory devices. They are used in laboratories in the chemical industry, schools and pharmacies to name a few. The devices are suitable for simultaneous mixing of substances which are placed on mounting plates. The built-in stirrer makes it possible to mix 5, 10 or 15 stirring positions at the same time and with equal intensity. The mixing intensity depends on the motor speed and the size of the magnetic stirring bar.

The model variants IKAMAG RT 5 power, IKAMAG RT 10 power and IKAMAG RT 15 power are also suitable for controlling temperature.

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers: an IKAMAG RO 5, 10 or 15 power and operating instructions bzw. an IKAMAG RT 5, 10 or 15 power and operating instructions.

Commissioning


Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

Function Stirring

The device is switched ON/OFF by means of the appliance switch **(A)** located on the right. The LED for "Power" **(B)** becomes green when the switch is activated.

The speed of the motor is set on the control knob "Mot" **(C)** with the scale from 1 to 10 **(D)**, which roughly corresponds to a range from 0 to 1200 rpm. The maximum speed of 1200 rpm is obtained when the control knob is turned all the way to the right in no-load operation.

If the viscosity of the medium increases during processing, then the speed drops slightly. If the viscosity decreases, e.g. due to the addition of low-


viscosity media, then the speed rises slightly. In either case the control knob can be adjusted to compensate for this speed change, provided the required torque does not exceed the maximum torque. Note that fluctuations in the mains voltage within the permissible tolerances will also result in slight variations in speed.

Function Heating

The device is switched ON/OFF by means of the appliance switch **(A)** located on the right. The LED for "Power" **(B)** becomes green when the switch is activated.

The temperature of the heater plate can be adjusted on the „Temp“ turn dial **(E)** with the corresponding scale from 1 to 10 **(F)**.

If energy is being directed to the heater plate, the „Temp“ LED **(G)** lights up yellow. The scale values 1 through 10 correspond to temperatures

on the surface from 0°C to 120°C.

Maintenance and cleaning

The IKAMAG RO 5 power, RO 10 power or RO 15 power and IKAMAG RT 5 power, RT 10 power or RT 15 power are maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Associated standards and regulations

Construction in accordance with the foolwing safety standards

EN 61 010-1 UL 3101-1 CAN/CSA C22.2 (1010-1)

Construction in accordance with the foolwing EMC standards

EN 50 082-1 EN 55 014 EN 60 555-2; -3

Associated EU guidelines

EMC-guidelines: 89/336/EWG

Machine guidelines: 73/023/EWG

Accessories

IKAFLOn stirrers, PTFE, 8 mm dia length 30mm

RSE Stirrer remover

Recommended IKA accessories

EH 4 basic suspension thermostat

RO 5.1 bath vessel

RO 10.1 bath vessel

RO 15.1 bath vessel

Technical data

Design voltage:		VAC	230±10%
	or	VAC	115±10%
Design frequency:		Hz	50/60
Motor:			ball-bearing brushless EC motor
Speed range:		rpm	0 to 1100 infinitely variable
Deviation between the - stirrer points		%	5
Perm. ambient temperature:		°C	+5 to +40
Perm. relative humidity:		%	80
Perm. ON duration:		%	100
Interference suppression:			to VDE 0871
Protection to DIN 40 050:			IP 42
Mounting surface		mm	120 x 440 (RO 5)
		mm	180 x 440 (RO 10)
		mm	270 x 440 (RO 15)
RT: Max surface temperature		°C	120
RT: Temperature constancy in the medium		K	± 2 depending on the container

List of spare parts

RO 5 / 10 / 15 power

RT 5 / 10 / 15 power

Item	Designation
1	Spacer
6	Deep-groove ball bearing
8	Protective cap
9	Plug
12	Clamping bush
13	Motor retainer
14	Potentiometer shaft
15	Rocker switch
16	BLP-Motor
17	EC - motor
18	Front plate
19	Hood
21	Bearing pins
22	Reflector
25	Miniature fuse
27	Locking ring
28	Pedestal
29	Driving belt
30	BLP Motor (small)
40	Cable clamp
41	Cable binder
42	Reference plate
52	Snap cable feed-through
53	Heater hood
2001	Magnet holder bearing pedestal (complete)
2002	Magnet holder motor (complete)

Max. stirring volume (water)		200 ml per stirring position in 250 ml Erlenmeyer flask, for other vessels
		2ltr (RO 5)
		4ltr (RO 10)
Stirring points		6ltr (RO 15)
		5 (RO 5)
		10 (RO 10)
Stirring point distance		15 (RO 15)
	mm	90 (RO 5)
	mm	90x90 (RO 10)
Max. power consumption:	mm	90x90 (RO 15)
	W	10 (RO 5)
	W	20 (RO 10)
Max. power consumption:	W	30 (RO 15)
	W	185 (RT 5)
	W	395 (RT 10)
Dimensions: (B x T x H)	W	600 (RT 15)
	mm	122 x 552 x 65 (RO 5)
	mm	182 x 552 x 65 (RO 10)
	mm	272 x 552 x 65 (RO 15)
	mm	138 x 552 x 65 (RT 5)
	mm	198 x 552 x 65 (RT 10)
Weight:	mm	288 x 552 x 65 (RT 15)
	kg	2,3 (RO 5)
	kg	3,2 (RO 10)
	kg	4,7 (RO 15)
	kg	3,0 (RT 5)
	kg	4,2 (RT 10)
	kg	6,0 (RT 15)

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