# IKA A10

# IKA M20

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*OPERATING INSTRUCTIONS* GB 7

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D

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## CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60 010-1; EN 50 081; EN 50 082; EN 55 014 and EN 60 555.

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
## Garantie

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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## Contens

	Page
Guarantee	2
Safety instructions	7
Correct use	7
Unpacking	7
Commissioning	7
Motor Circuit Breaker	8
Description	8
Grinding	8
Grinding while cooling	9
Replacing the grinding tool	9
Maintenance and cleaning	9
Associated standards and regulations	10
Accessories	10
Technical data	10
List of spare parts A10	16
Spare parts diagramm A10	17
List of spare parts M20	18
Spare parts diagramm M20	19

## Safety instructions

The A10 or M20 devices may not be used in areas subject to explosion hazards.

Protective glasses must be worn when performing experimental cuts. In addition, suitable protective gloves must be worn when adding coolants.

Experiments involving comminution must be performed behind a protective wall (Plexiglas®-wall, exhaust hood or the like) in order to prevent injury caused by the sudden occurrence of excess pressure (possible bursting of the cover).

Please comply with all of the relevant safety regulations.

## Correct use

IKA analysis grinders A10 and IKA universal grinder M20 are suitable for carrying out low-loss grinding of hard and brittle substances as of grain size 6 - 7mm. The grinders operate without creating dust or lost. Beaters, blades and grinding are replaceable. The useful volume of the grinding chamber is 50ml for the A10 and 250ml for the M20 grinder.

Hard and brittle ground material is crushed by means of beaters or blades rotating at high speed.

Oil substances or those that are given to smearing and which are not capable of being ground at room temperature have to be cooled for grinding (e.g. embrittled with liquid nitrogen).

## Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The A10 delivery quantity includes: one A10 grinder, one A14 stainless steel beater, one A18 grinding chamber reduction element, one set of tools and operating instructions.

The M20 delivery quantity includes: one M20 grinder, one M21 stainless steel beater, one set of tools and operating instructions.

## Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these proce-

dures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

**NOTE:** The grinder is designed for short term operation. Maximum operating time **A10** - 5 minutes and **M20** - 7 minutes.

After that, the grinder must cool off for approximately 10 minutes. This condition does not apply if the grinders is connected to a cooling unit.

## Motor Circuit Breaker

**A 10** The device is electronically protected against overload and excessive temperatures. A error yellow lamp lights on the front panel up if the fuse is triggered and the motor comes to a standstill.

**M 20** The grinder is protected against overload by means of an overcurrent protective device. In the event of overload the overcurrent protective device interrupts the current supply. Restarting is carried out by pressing-in the black button at the rear of the grinder. The motor is switched on via push-button at the front.

## Description

**A 10** The device consists of a drive unit with grinding chamber and a cover. The units are connected to each other by means of quick-release locks. When assembling the two units make certain that the safety pin in the drive unit is pushed down by the pressure pin on the cover. Otherwise, a tight and reliable connection is not ensured, the yellow error lamp will light up and the grinder will not start. The device is started by pressing the push-button at the front.

**M 20** The device consists of three basic units (drive unit, middle part and cover) which are connected to each other by means of quick-release locks. When assembling the three units make certain that the hose coupling for cooling the middle part is located at the rear of the device. When putting on the cover make certain that the safety pin in the middle part is pushed down by the pressure pin on the cover. Otherwise, a tight and reliable connection is not ensured and the grinder will not start. The device is started by pressing the push-button at the front.

## Grinding

Fineness is mainly dependent on how long the grinding operation is performed.

**A 10** Dry, brittle materials are easily ground, but may not be larger than the size of a pea as there may be danger of damage when carbide metal beaters are used. The A14 stainless steel beater is assembled from rust- and acid-proof steel (material no. 1.4034). The A15 carbide metal beater made of carbide metal (86% tungsten carbide and 14% cobalt) is available if particularly wear-resistant material is required. The A17 star-type blade (material no. 1.4310) is available for specifically light ground materials such as pharmaceuticals, tablets, cardboard, felt, polystyrene, plastic foams (Brittle) and the like. The A18 grinding chamber reduction element (included in the delivery quantity) is inserted into the grinding chamber for comminution of small sample quantities (1-5ml) and specifically lighter products.

**M 20** If materials are to be ground that are smaller than a certain grain size then it should be sieved and the residue in the sieve is put back into the grinder. For sharp-edged material which is thus very abrasive, it is recommended that the entire quantity is left in the grinder in order to additionally exploit the properties of the material.

## Grinding while cooling

Plastic, soft or smeary materials have to be cooled if they are to be ground.

There are two hose couplings located at the rear of the grinder in order to connect the cooling brine. Cooling with tap water is usually sufficient for substances which become hot as a result of the grinding operation and are thus given to smearing or those that produce oil.

The grinding chamber may only be filled halfway especially when drugs containing oils and similar materials are ground.

If more refrigeration is required then the ground matter should be cooled prior to being put into the grinding chamber.

The following procedure shall apply if liquid nitrogen is used in order to embrittle the material to be ground:

### **Safety information:**

Suitable protective gloves must be worn when filling and transferring liquid nitrogen!

Protective glasses must be worn!

- The sample to be ground is put into a suitable vessel (e.g. evaporating dish) and covered with enough liquid nitrogen until it is completely covered.

When using liquid nitrogen make certain that the sample is not capable of reacting with oxygen.

- Cool grinding chamber using liquid nitrogen.
- Put the deep-cooled sample into the grinder and grind approximately 15 seconds.
- Then open the grinder and add liquid nitrogen for only as long as it takes for the sample to become deep-cooled again (embrittled).

### **Safety information:**

Make certain that the liquid nitrogen has evaporated before beginning the grinding operation.

NOTE: Danger of overpressure!

- Repeat the grinding operation for approximately 15 seconds.

- Repeat the sequence of cooling and grinding until desired fineness has been achieved.

- If the required fineness has not been achieved after 5 minutes, then a further reduction in grain size is unlikely.

The respective grinding chambers are sealed by means of labyrinth sealings. Though they are dust-proof, they are not leak-proof. This should be taken into account for substances where liquids may escape.

## Replacing the grinding tools

The grinding tool is fixed by means of the retaining tube and the cap nut is unscrewed by using the SW8 or 10 spanner. The grinding tool, which rests on a double-edge track, may then be pulled out.

## Maintenance and cleaning

The IKA A10 and M20 are maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

**When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.**

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical equipment may not be put into the cleansing agent for cleaning.

The operator shall consult the manufacturer before using a method of cleaning or decontamination other than the one which is recommended by the manufacturer in order to make certain that the method which the operator intends to employ will not damage the equipment.

## Associated standards and regulations

### Associated EU guidelines

EMC guideline: 89/336/EWG, 93/31/EWG

Low voltage guideline: 73/023/EWG

Machine-guideline: 89/392/EWG, 93/44/EWG, 93/68/EWG

### Construction in accordance with the following safety standards

EN 61 010-1

EN 292-1, -2

UL 3101-1

EN 414

CAN/CSA C22,2 (1010-1)

### Construction in accordance with the following EMV standards

EN 55 014-1

EN 50 081

EN 50 082-1

## Accessories

### A 10

**A 18** grinding chamber reduction element

**A 11** cover

**A 14** spare beater, stainless steel

**A 15** beater, carbide metal

**A 17** star-type blade, stainless steel

### M 20

**M 20.1** spar grinding chamber with cover

**M 21** spare beater, stainless steel

**M 22** beater, stainless steel

**M 23** star-type blade, stainless steel

## Technical data

Design voltage:	<b>VAC</b>	230±10%	
	or	<b>VAC</b>	115±10%
Design frequency:	<b>Hz</b>	50/60	
Permissible ambient temperature:	<b>°C</b>	+5 bis +40	
Permissible relative humidity:	<b>%</b>	80	
Degree of protection			
as in DIN 40 050:		IP 21	
Rated speed:	<b>1/min</b>	20 000	
Motor:		universal motor	

### **A 10**

Maximum motor input:	<b>W</b>	180
Grinding chamber useful volume:	<b>ml</b>	50
Permissible operating time		
(without cooling):	<b>min</b>	5
Interference suppression:		as in VDE 0871
Dimensions:	<b>mm</b>	120 x 105 x 225
Weight:	<b>kg</b>	2,2

### **M 20**

Maximum motor input:	<b>W</b>	550
Grinding chamber useful volume:	<b>ml</b>	250
Permissible operating time		
(without cooling):	<b>min</b>	7
Interference suppression:		as in VDE 0871
Dimensions:	<b>mm</b>	170 x 170 x 350
Weight:	<b>kg</b>	6,6

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## List of spare parts

### A 10

Item	Designation
7	cup nut
9	grooved pin
13	distance ring
14	O-ring
15	O-ring
16	Nilos-ring
17	Nilos-ring
19	ball bearing compensating disc
20	fan blade
21	cap
22	protective rubber sleeve
25	housing
26	flange
27	felt ring
28	hub
29	threaded bushing
31	holder
32	foundation
34	front plate
39	spar cutter A14
40	grinding bowl
43	spring
45	eccentric lever
48	commutator motor
50	printed circuit board
51	cover A11
52	connection line
2001	grooved ball bearing
2002	carbon brush
3001	push-button
3002	push-button

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




---

LABORATORY TECHNOLOGY

ANALYSING TECHNOLOGY  
PROCESSING EQUIPMENT

**IKA® WORKS, (Asia) Sdn Bhd**

---

LABORATORY TECHNOLOGY

ANALYSING TECHNOLOGY  
PROCESSING EQUIPMENT

Europa - Afrika

---

**IKA®** LABORTECHNIK  
JANKE & KUNKEL-STR. 10  
79 219 STAUFEN  
GERMANY  
TEL. 07633/831-0  
FAX 07633/831-98  
E-mail: [sales@ika.de](mailto:sales@ika.de)

America

---

**IKA®** WORKS, INC.  
2635 NORTH CHASE PKWY, SE  
WILMINGTON, NC 28405-7419  
TEL. 800/733-3037  
TEL. 910/452-7059  
FAX 910/452-7693  
E-mail: [ikaworks@isaac.net](mailto:ikaworks@isaac.net)

Asien - Australien

---

**IKA® Works (Asia) Sdn Bhd**  
(Company No. 340448-K)  
No. 3, Lot PT 3403,  
Batu 22<sup>3/4</sup> Lebuhraya,  
Kuala Lumpur-Serendah,  
48200 Serendah  
Selangor, Malaysia  
TEL : (603) 601 3122  
FAX : (603) 601 3126  
E-mail: [ika@tm.net.my](mailto:ika@tm.net.my)