IKA® RW 47 D

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
OPERATING INSTRUCTIONS GB 11

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Contens

	Page
Safety instructions	11
Correct use	12
Unpacking	12
Useful facts	13
Drive	13
Motor circuit breaker	13
Speed display	14
Gear speed change	14
Commisioning	14
Switching on the machine	15
Output shaft	15
Maintenace and cleaning	15
Technical data	16
Accessories	17
Permitted IKA stirrer tools	17
Associated standards and regulations	17
Guarantee	17
Verdrahtungsplan	77
Spare parts list	78
Spare parts diagram	79

Safety instructions

Never operate the equipment with the stirrer tools rotating freely. Ensure that parts of the body, hair or items of clothing cannot be trapped by the rotating parts. Pay attention when setting the speed to any imbalance of the stirrer tools and possible spraying of the medium to be stirred. Use stirring shaft protective equipment and a safety cutout for height adjustable stands. The equipment is not suitable for manual operation. Please follow the relevant safety instructions and guideli-

nes, and occupational health and safety regulations for use in the laboratory. When working with the stirring unit, the user must select and wear his personal protective equipment according to the mixing hazard category. Defective or inappropriate protective equipment can expose the user to the risk of splashing or getting caught in the stirring element.

It is important to note that the surfaces of the motor (cooling fins) and certain parts of the roller bearing can get very warm during operation.

Avoid knocks and impacts on the lower end of the shaft and the chuck gear teeth. Even minor, invisible damage can lead to imbalance and uneven shaft action. Careful handling guarantees safe working and the long life of the equipment.

Imbalance of the output shaft, the chuck and in particular the stirring tools can lead to uncontrolled resonant vibrational behaviour of the equipment and the whole assembly. Glass apparatus and stirrer containers can be damaged or shattered by this. It can cause injury to the operator, as can the rotating stirring tool. If vibration of the equipment is noticed, the speed must be reduced immediately in all cases until no more vibration occurs. In this case exchange the stirring tool for one without imbalance or remedy the cause of the imbalance.

Operation with a freely rotating stirring element without a stirring vessel is dangerous. The stirring elements may only be inserted into the hollow shaft to such a depth that they do not abrade the top edge of the casing. Ensure that the stirring element is securely clamped in the chuck.

The unit must be switched off before changing to another gear speed, otherwise the toothed gears in the reduction stage will be damaged.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening. Live parts inside the machine may still be live for some time after unplugging from the mains.

The high torque developed by the RW 28 basic requires particular care in the choice of stand, cross sleeve and anti-rotation element for the stirring bowl.

NOTE! Covering or parts that are capable of being removed from the unit without accessory equipment have to be reattached to the unit for safe operation in order to prevent, for example, the ingress of fluids, foreign matter, etc.

The agitated vessels used for stirring have to be secured.

If the stirring unit is operated in such a way that the appliance's EMERGENCY OFF switch is not accessible, you must add another EMERGENCY OFF switch to the work area that can be easily reached.

Additional hazards to the user may occur if inflammable materials are used during the stirring operation or, if glass equipment is used, the energy transmitted from the stirring device results in glass breakage.

Make certain that the unit is set at the lowest speed before commissioning; otherwise, the unit will begin running at the speed last set. Stirring at high speeds can cause the medium to splash out of the mixing vessel. There may be dangerous electrostatic activity between the medium and the output shaft.

The appliance does **NOT** start up again automatically following a cut in the power supply or mechanical failure during a stirring procedure.

Stirring unit RW 47 D may not be used in explosive atmospheres, for mixing dangerous substances or under water.

Do not touch rotating parts during operation!

Only use IKA accessories.

Accessories may only be assembled once the plug has been disconnected.

Use stirring shaft protective equipment and a safety cutout for height adjustable stands.

Correct use

Laboratory stirrer RW 47D can be used to stir and mix liquids of low to medium viscosity with various stirring tools. It is designed for use in technical colleges and for small-scale production. The appliance must be carefully secured to a stand with four screws for proper use (see **fig. 2**). You must also ensure that the mounting equipment (e.g. floor stand, telescope stand) is set up in such a way that it is not liable to topple and does not start to move during the stirring procedure.

Unpacking

Please unpack the machine carefully and inspect for damage. It is important that any transit damage should be noted at the time of unpacking. In certain circumstances it may be necessary to investigate immediately (post, rail or freight forwarder).

The guarantee covers:

An RW 47 D stirring unit as per order, a chuck key, allen key, allen screws and operating instructions.

Useful facts

By buying this appliance, you have chosen a high quality product. The appliance has a sturdy construction and is simple to operate, thus ensuring ease of use. You can set speeds of 57 to 1300 1/min with the 50Hz (mains frequency) version and speeds from 69 to 1560 1/min with the 60Hz (mains frequency) version. The unit heats up during use. The ample cooling surfaces on the motor serve to distribute and radiate the heat as evenly as possible.

You can also use high-sided mixing vessels together with IKA stands without any problem.

The unit comes without a connection cable because the user's installations mean that it is not possible to issue a standard-length cable due to the different sizes of premises. Experts can determine the pin assignment from the wiring plan.

The precision shaft is shaped as a hollow shaft and accommodates stirrer shanks of up to 13 mm in diameter. It is important to note here that these may only be inserted to just before the top edge of the casing. You must also ensure that the stirrer shank is securely clamped in the chuck.

Drive

The method of operation of the motor is characteristic of the operating behavior of the RW 28 basic units. The motor can be infinitely adjusted via a friction wheel drive unit but is operated in one position only. The power output, rotational speed and torque of the motor are regarded as constant and are optimized for this operating position. After the friction wheel drive the motor power output is transmitted to the off-drive shaft of the stirrer via either the first or second stage of the subsequent manually adjustable toothed gear train. Neglecting all losses, the power at the off-drive shaft during

power transmission (conversion) is always as great as the power output of the motor. The gear train merely provides rpm and torque conversion, which can be described in purely theoretical terms by the characteristic curve plotted in Fig. 1.

Friction wheel wear is kept low by a helical clutch matched to the gear train. The down-pressure required at the friction wheel is adjusted by the helical clutch according to the moment acting on the stirrer shaft. A low torque results in a low down-pressure; a high torque results in a high down-pressure.

Motor circuit breaker

The ball race three-phase motor is maintenance-free. Its attached switching unit features a special safety device for operation in the case of overcurrent and undervoltage (undervoltage release) to permanently disconnect the motor and avoid thermal damage. The appliance can only be started up again by activating the ON/OFF switch (**fig. 2**), when power permitted for the motor is supplied.

You can find the basic connection and rating for the motor on its rating plate.

The safety switch SI 400 can also be connected in the motor switching unit. When used in combination with an IKA stand, this allows the user to switch on the appliance only at a defined height in the stand's working range.

If the stirring unit is operated in such a way that the appliance's EMERGENCY OFF switch ((B)-**fig. 2**) is not accessible, you must add another EMERGENCY OFF switch to the work area so that it can be easily reached.

Speed display

The speed on the RW 47 D agitator is adjusted by means of the rotary button (C).

You can read the set speed from the marking on the top of the casing on the label attached to the motor support. (The values are only valid for low agitator shaft lows).

Gear speed change


You can select one of the three gear positions available (0, I or II) using the toggle button (D) on the front of the appliance. You may only switch to a different cog wheel gear level when the appliance is switched off (motor stopped). Otherwise this may cause damage to the gear cogs. The arrow on the front plate indicates the direction of the chuck.

Commissioning

Laboratory stirrer RW 47D must be secured to a stable stand with four screws (e.g. R474 or R472) and wired according to the wiring plan for proper use. The stirring bowl must always be well secured for safety reasons. You must also ensure that the mounting equipment (e.g. floor stand, telescope stand) is set up in such a way that it is not liable to topple and does not start to move during the stirring procedure. Accessories must be assembled and wired according to the assembly instructions.

If the stirring unit is operated in such a way that the appliance's EMERGENCY OFF switch ((B)-**fig. 2**) is not accessible, you must add another EMERGENCY OFF switch to the work

area so that it can be easily reached. Points of particular danger (such as crushing, rotating parts, clamping screws for stand legs and rollers etc.) are marked in **fig. 3** with an exclamation mark. Please pay particular attention to these danger points during operation!

Switching on the machine

Check whether the voltage given on the type plate corresponds to the available mains voltage. The socket used must be earthed (fitted with earth contact). If these conditions have been met, the machine is ready to operate when the mains plug is plugged in. If these conditions are not met, safe operation is not guaranteed and the machine could be damaged.

When working with safety switch SI 400, the appliance only starts up once it has been set to the designated working height and the switch contact activates the safety switch.

During commissioning of the machine the output shaft starts to run at the last speed set. Therefore check the setting of the control knob (speed). Also ensure that the speed set is suitable for the test texture selected. If you are unsure, set the motor slide to the lowest speed (front motor slide stop) using the side knob and the gear to the lowest level (**fig. 2**).

The appliance will start up when you press the ON switch on the motor (**fig. 2**).

If the unit has not been used for some time, a knocking noise will be heard when switching on, which is due to the preload on the friction lining of the friction wheel drive. This has no detrimental effect on the function of the unit, and the knocking noise will disappear after a short running-in period.

Output shaft

The chuck (E) - **fig. 2**) and off-drive shaft allow stirring elements with shank diameters of 3mm to 16mm to be clamped. Shank diameters of up to 13mm can be inserted through the hollow shaft to just before the top edge of the inside casing. You must ensure that the stirring shank is securely clamped in the chuck.

Maintenance and cleaning

The stirrer is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

When cleaning the unit, take particular care not to bring substances containing oil or grease into contact with the running surfaces of the cone pulley and the friction wheel. This would substantially reduce the frictional coefficient resulting from the pairing of the friction wheel and cone pulley materials, and would impair power transmission in the unit.

Make certain that the hexagon head nut (part no. 153 of spare parts diagram) is sufficiently tightened in order to prevent the shifting of the motor plate and thus the speed as well.

Technical data

Speed range

50Hz stage I	min ⁻¹	57 - 275
50Hz stage II	min ⁻¹	275 - 1300
60Hz stage I	min ⁻¹	69 - 330
60Hz stage II	min ⁻¹	330 - 1560
Max. torque stirrer shaft (Overload)	Ncm	4642 (with 57 1/min)
Max. torque stirrer shaft	Ncm	3000 (with 100 1/min)
Max. torque stirrer shaft	Ncm	285 (with 1000 1/min)
Permitted on-time:	%	100
Speed adjustment:	-	mechanical rotary button
Speed display:	-	on the side speed label

Motor data

Design voltage:	VAC	3x400 Y	3x230 Δ	3x400 Y	3x230 Δ
Design frequency:	Hz	50	50	60	60
Motor speed:	min ⁻¹	2880	2850	3425	3425
Input power:	W	554	548	489	490
Output power:	W	370	370	370	370
Maximum output at stirrer shaft:	W	300			
Overall efficiency:	%	54			
Protection type to DIN EN 60529:	-	IP 44			
Excess voltage category:	-	II			
Protection class:	-	I (protective earth)			
Contamination level:	-	2			
Protection at overloaded:	-	overcurrent and undervoltage switch on the three-phase motor			
Perm. ambient temperature:	°C	+5 to+40			
Perm. ambient humidity: (rel.)	%	80			
Operating position:	-	on stand, clamping chuck pointing down			
Drive:	-	Rib cooled three-phase motor with friction wheel drive and subsequent 2-stage toothed gear train			
Maximum stirring quantity - water:	ltr	200			
for high-viscosity (HV):	mPas	10000-100000			
Clamping chuck - clamping range:	mm	3 - 16			
Hollow shaft internal diameter:	mm	13 (to short before the internal housing upper edge)			
Stand attachment:	-	Flange on housing lower part			
Housing:	-	Al, hermoplastic			
Dimension without extension arm: (WxDxH)	mm	145x445x340			
Weight	kg	15			
Operation at a terrestrial altitude	m	max. 2000 above sea level			

Accessories

R 472	Floor stand
R 474	Telescopic stand
R 302	Stirring shaft protection
SI 400	Safety switch
SI 472	Holding fixture
SI 474	Holding fixture

Permitted IKA stirrer tools

Max. Speed 1/min

R 2305	Propeller stirrer 4fl	1 300
R 2302	Propeller stirrer 4fl	600
R 1385	Propeller stirrer 3fl	800
R 1388	Propeller stirrer 3fl	400
R 1313	Turbinenrührer	800
R 1300	Dissolver stirrer	2 000
R 1301	Dissolver stirrer	800
R 1302	Dissolver stirrer	1 000
R 1303	Dissolver stirrer	2 000
R 1376	Surface stirrer	800
R 2311	Surface stirrer	600
R 1333	Anchor stirrer	800
R 1393	Worm stirrer	1 000

Associated standards and regulations

Associated EU guidelines

EMC-guidelines:	89/336/EG	92/31/EG
Low voltage guidelines:	73/023/EG	
Machine guidelines:	98/37/EG	

EU standards

IEC 61010-1	DIN EN 61 010-1	EN 60 204-1
EN 292-1, -2	EN 414	

Int. standards

UL 3101.1	CAN/CSA C22.2 (1010.1).
IEC 61 326 –1	IEC 61 326 –1

Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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

spare parts

Item	Designation	Item	Designation
3	chuck	124	rotary control knob
6	strip	152	cone pulley
15	warning "attention"	153	threaded pin
21	feather	154	rack pinion
23	grooved ball bearing	155	speed label scala
24	grooved ball bearing	156	distance bushing
31	cylindrical roller bearing	157	pressure element
32	deep-groove ball thrust b.	158	pressure spring
33	needle bearing	159	marking label
34	steel ball	162	conveyor plate
35	radial shaft-sealing ring	167	motor support
43	Nilos ring	174	housing-upper part
48	bearing plate	2001	friction wheel, compl.
49	off-drive shaft	2002	gear unit, compl.
50	ageing washer	2003	motor, lacquered
51	thrust washer	5001	grooved ball bearing
52	drive gear	5002	grooved ball bearing
53	coupling upper part	5003	plastic hood
54	coupling bottom part	5004	fan blade
55	pressure spring	5011	contactor modell D
56	pressure spring	5012	housing ground
57	reverse label		(motor starter)
58	counter gear	5013	housing upper part
59	driven gear		(motor starter)
60	toggle	5014	callipers I ON green
66	housing bottom part	5015	callipers 0 OUT/RESET
72	switch shaft		red
79	friction wheel		
97	packing gland bolting		
98	nut PG-bolting		
101	three-phase motor lacquered		
102	motor starter incl. housing		
103	overload relay		

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.

2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd

No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.

293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.

Estrada do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-9600
FAX +55 21 2435-9601
E-mail: brasil@ika.net

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou

173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn