



MF 10.2 Impact grinding head

For crushing brittle, hard materials such as minerals, building materials up to Mohs hardness 6. Before being discharged, the ground material passes through a sieve. This sieve is interchangeable and available in different hole sizes (not incl. with delivery). The ground material can then be collected using an NS 29 standard ground vessel.

Technical Data

Process type	inline
Circumferential speed max. [m/s]	31.4
Feed hardness max. [Mohs]	6
Feed grain size max. [mm]	10
Material milling chamber	stainless steel 1.4308
Milling chamber, can be cooled with water	no
Mill feed can be cooled in milling chamber with dry ice	yes
Mill feed can be cooled in milling chamber with liquid nitrogen	no
Dimensions (W x H x D) [mm]	320 x 450 x 300
Weight [kg]	2.7
Ident. No.	2871000