



PROMOTION



Rotary Evaporators

"all in one" Packages

RV 10



RV 10 control V "all in one" Package	
Model	Ident. No.
RV 10 control V (includes Heating bath HB 10 control and set of glassware vertical)	8022700
MPC 105 T Vacuum pump	4067700
KV 600 digital Chiller	3410500
RV 10.4002 Magnetic valve	3880300

RV 10 control V "all in one" Package price
 Package list price:
~~USD \$ 10,668.00~~

Promotion price USD \$ 8,934.00

-16%



RV 10 digital V "all in one" Package	
Model	Ident. No.
RV 10 digital V (includes Heating bath HB 10 digital and set of glassware vertical)	8022500
LVS 105 T 10-ef Vacuum pump	4067800
KV 600 digital Chiller	3410500

RV 10 digital V "all in one" Package price
 Package list price:
~~USD \$ 12,260.00~~

Promotion price USD \$ 9,873.00

-19%



RV 10 control V

- With integrated vacuum controller
- Fully automated rotary evaporator
- Complete with volume dependent distillation programs
- Integrated upgradable solvent library
- Color graphic display for safe and comfortable operation
- Display of distillation curves
- Heating bath safety management
- USB interface
- Complete with glassware set

RV 10 digital V

- Digital water / oil heating bath with integrated carrying handles
- Digital temperature display
- Infrared interface for data transfer
- Complete with glassware set

MPC 105 T Vacuum pump

- For dry, oil-free applications in the low vacuum range
- Compact and space saving design
- Quiet running and long lifetime
- Chemically resistant diaphragm pumps (MPC), are resistant to aggressive solvents & acidic vapors. The diaphragms and gas contacting parts consist of PTFE & PTFE compounds. The pump & connection heads are carbon fiber reinforced with electrical conductivity.

LVS 105 T 10-ef Vacuum pump

- Compact, dry-running, chemical resistant pump systems, the ideal solution for many applications in chemical laboratories and research
- Fully automatic pump system with built-in solvent library
- Graphic display
- Close to 100% solvent recovery
- Equipped with a speed regulated vacuum controller and a chemical resistant diaphragm pump. The pumping speed can be precisely controlled, e.g. for gentle distillations

KV 600 digital Chiller

- KV 600 digital is an actively cooled chiller with an air-cooled refrigeration system
- Featuring a user-friendly microprocessor controller with large temperature display
- The temperature consistency is 1 K



Designed to work perfectly

IKA® -Works (Asia) Sdn Bhd
 No. 17 & 19, Jalan PJU 3/50
 Sunway Damansara Technology Park
 47810 Petaling Jaya Selangor, Malaysia

Tel. +60 3 7804-3322
 Fax +60 3 7804-8940
 sales@ika.com.my

Price are on ex-works Guangzhou, China basis and subject to additional importation costs and local taxes. We reserve the right to replace these prices, by new, subsequent ones as a response to changing circumstances, this right also applies if the printed price is a blatantly obvious printing error. Prices valid until 31 October 2010. Subject to technical change.